

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

SSCP-013

Pro-PICC® Power Injectable Peripherally Inserted Central Catheter

IMPORTANT INFORMATION

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device.

This SSCP is not intended to replace the Instructions for Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients.

Applicable Documents	
Document Type	Document Title / Number
DHF	10004, 11010-A1, 11010-A3, 11010, 11011-A3, 11011, 11012-A1, 11012
'MDR Documentation' File Number	MDR-013

Revision History					
Revision	Date	CR#	Author	Description of Changes	Validated
1	25APR2022	26921	RS	Implementation of SSCP	<input type="checkbox"/> Yes, this version was validated by the Notified Body in the following language: English <input type="checkbox"/> No, this version was not validated by the Notified Body as this is a Class IIa or IIb implantable device

Revision History					
Revision	Date	CR#	Author	Description of Changes	Validated
2	17JUN2022	27027	RS	Scheduled Update	<input type="checkbox"/> Yes, this version was validated by the Notified Body in the following language: English <input type="checkbox"/> No, this version was not validated by the Notified Body as this is a Class IIa or IIb implantable device
3	23NOV2022	27509	GM	Scheduled Update; updated SSCP in accordance with CER-013_C and QA-CL-200-1 Version 3.00 Template. Acronym table was added in Section 7 of the Patient Section.	<input checked="" type="checkbox"/> Yes, this version was validated by the Notified Body in the following language: English <input type="checkbox"/> No, this version was not validated by the Notified Body as this is a Class IIa or IIb implantable device
4	20OCT2023	28545	GM	Update in accordance with CER-013_D	<input type="checkbox"/> Yes, this version was validated by the Notified Body in the following language: English <input type="checkbox"/> No, this version was not validated by the Notified Body as this is a Class IIa or IIb implantable device
5	24OCT2024	29499	GM	Update in accordance with CER-013_E	<input type="checkbox"/> Yes, this version was validated by the Notified Body in the following language: English <input type="checkbox"/> No, this version was not validated by the Notified Body as this is a Class IIa or IIb implantable device

USERS / HEALTHCARE PROFESSIONALS

The following information is intended for users/healthcare professionals. Following this information there is a summary intended for patients.

1. Device identification and general information

Device trade name(s)	Pro-PICC® Power Injectable Peripherally Inserted Central Catheter
Manufacturer name and address	Medical Components, Inc. 1499 Delp Drive Harleysville, PA 19438 USA
Manufacturer single registration number (SRN)	US-MF-000008230
Basic UDI-DI	00884908286NP
Medical device nomenclature description / text	C010201 – Central I.V. Catheters, Peripheral Access
Class of device	III
Date first CE certificate was issued for this device	Pro-PICC®- October 2007 Pro-PICC® Valved - May 2013 Jet-PICC - July 2009 PFM-PICC - December 2015
Authorized representative name and SRN	Gerhard Frömel European Regulatory Expert Medical Product Service GmbH (MPS) Borngasse 20 35619 Braunfels, Germany SRN: DE-AR-000005009
Notified Body name and single identification number	BSI Group The Netherlands B.V. NB2797

The devices in scope of this document are all peripherally inserted central catheter (PICC) sets. The device part numbers are organized into variant categories. These devices are distributed as procedure trays, in various configurations inclusive of accessories and adjunctive devices (see section “Accessories intended for use in combination with the Device”).

Variant Devices:

Variant Description	Part Number(s)	Explanation of Multiple Part Numbers
3F x 55cm Single Lumen Pro-PICC®	10467-855-800 10467-855-801	No significant clinical, biological, or technical difference (only difference is branding)
4F x 55cm Single Lumen Pro-PICC®	10602-855-800 10602-855-801	No significant clinical, biological, or technical difference (only difference is branding)
4F x 55cm Single Lumen Valved Pro-PICC®	10643-855-801	N/A
5F x 55cm Double Lumen Pro-PICC®	10561-855-800 10561-855-801	No significant clinical, biological, or technical difference (only difference is branding)
5F x 55cm Double Lumen Valved Pro-PICC®	10645-855-801	N/A
5F x 60cm Single Lumen Pro-PICC®	10556-860-800 10556-860-801	No significant clinical, biological, or technical difference (only difference is branding)
5F x 60cm Single Lumen Valved Pro-PICC®	10644-860-801	N/A
6F x 60cm Double Lumen Pro-PICC®	10563-860-800 10563-860-801	No significant clinical, biological, or technical difference (only difference is branding)
6F x 60cm Triple Lumen Pro-PICC®	10568-860-800 10568-860-801	No significant clinical, biological, or technical difference (only difference is branding)
6F x 60cm Triple Lumen Valved Pro-PICC®	10646-860-801	N/A

Procedure Trays:

Catalog Code	Part Number	Description
JSACT5D	10561-855-800	5F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
JSACT5DL	10561-855-800	5F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
JSACT6D	10563-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
JSACT6DL	10563-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
PFMCT5DLWS	10561-855-800	5F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT5DS	10561-855-800	5F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP52024	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP52028	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE

Catalog Code	Part Number	Description
MRCTP62024	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP62028	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MR17035201	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR17035202	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
MR17035205	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN NURSING SET
MR17036201	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR17036202	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
MR17036205	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN NURSING SET
JSACT4S	10602-855-800	4F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
JSACT4SL	10602-855-800	4F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
JSACT5S	10556-860-800	5F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
JSACT5SL	10556-860-800	5F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
PFMCT3SS	10467-855-800	3F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
PFMCT4SLWS	10602-855-800	4F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT4SS	10602-855-800	4F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
PFMCT5SLWS	10556-860-800	5F X 60CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT5SS	10556-860-800	5F X 60CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP31024	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP41024	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP41028	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MRCTP51024	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP51028	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MR17033101	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET

Catalog Code	Part Number	Description
MR17033102	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17033105	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
MR17034101	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR17034102	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17034105	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
MR17035101	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR17035102	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17035105	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
JSACT6T	10568-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET
JSACT6TL	10568-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN LONG WIRE SET
MRCTP63024	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN FULL CUT DOWN SET
MR17036301	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET
MR17036302	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN LONG WIRE SET
MR17036305	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN NURSING SET
MR82034101	10643-855-801	4F X 55CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR82035101	10644-860-801	5F X 60CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR82035201	10645-855-801	5F X 55CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR82036301	10646-860-801	6F X 60CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET

Configurations of Procedure Trays:

Configuration Type	Kit Components
Pro-PICC® Basic Set	(1) Catheter w/ Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (3F Sets) 1.1mm ID x 10cm (3.5F) Peelable Introducer, (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement

Configuration Type	Kit Components
	Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Pro-PICC® Longwire Set	(1) Catheter, (1) Peelable Introducer: (3F Sets) 1.1mm ID x 10cm (3.5F) Peelable Introducer, (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 130cm (.018) Coated Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Pro-PICC® Nursing Set	(1) Catheter w/ Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Tourniquet, (1) Peelable Introducer: (3F Sets) 1.1mm ID x 7cm (3.5F) Peelable Introducer, (4F Sets) 1.5mm ID x 7cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 7cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 7cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Safety Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 45cm (.018) Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Pro-PICC® Full Cut Down Set	(1) Catheter w/Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (3F Sets) 1.1mm ID x 7cm (3.5F) Peelable Introducer, (4F Sets) 1.5mm ID x 7cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 7cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 7cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Safety Needle w/ Echo Tip, (1) 0.47mm x 45cm (.018) Guidewire Floppy Straight Tip, (1) Scalpel, (1) Securement Device, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Pro-PICC® Full Cut Down Set with 70cm Guidewire	(1) Catheter w/ Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip, (1) Scalpel, (1) Securement

Configuration Type	Kit Components
	Device, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Pro-PICC® Valved Basic Set	(1) Catheter, (1) Peelable Introducer: (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient Information Packet, (1) Patient ID Card
Jet-PICC CT Basic Set	(1) Catheter w/ Stylet, (1) 0.76mm (0.030) I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
Jet-PICC CT Longwire Set	(1) Catheter, (1) Peelable Introducer: (4F Sets) 1.5mm ID x 10cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 10cm (5.5F) Peelable Introducer, (6F Sets) 2.0mm ID x 10cm (6.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 10cc Syringe, (1) 0.47mm x 130cm (.018) Coated Guidewire Floppy Straight Tip, (1 2 3) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet
PFM-PICC Full Cut Down Set	(1) Catheter w/ Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (3F Sets) 1.1mm ID x 7cm (3.5F) Peelable Introducer, (4F Sets) 1.5mm ID x 7cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 7cm (5.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Safety Needle w/ Echo Tip, (1) 0.47mm x 45cm (.018) Coated Guidewire Floppy Straight Tip, (1 2) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet

Configuration Type	Kit Components
PFM-PICC Full Cut Down Set with 70cm Guidewire	(1) Catheter w/ Stylet, (1) 0.76mm (0.030") I.D. Adaptor w/ Sideport, (1) Peelable Introducer: (4F Sets) 1.5mm ID x 7cm (4.5F) Peelable Introducer, (5F Sets) 1.8mm ID x 7cm (5.5F) Peelable Introducer, (1) 0.9mm OD x 0.5mm ID x 70mm (21GA) Needle w/ Echo Tip, (1) 0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip, (1 2) Needleless Connector(s), (1) Securement Device, (1) Scalpel, (1) Tape Measure, (1) Patient ID Card, (1) Patient Information Packet

2. Intended use of the device

Intended purpose	<p>Pro-PICC®, Jet-PICC, and PFM-PICC Pro-PICC®/Jet-PICC/PFM-PICC Power Injectable Peripherally Inserted Central Catheters are intended for use in adult and pediatric patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician. The device is intended to be used under the regular review and assessment of qualified health professionals. This catheter is for Single Use Only.</p> <p>Pro-PICC® Valved Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheters are intended for use in adult patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician. The device is intended to be used under the regular review and assessment of qualified health professionals. This catheter is for Single Use Only.</p>
Indication(s)	The Pro-PICC®/Pro-PICC® Valved/PFM-PICC/Jet-PICC Power Injectable Peripherally Inserted Central Catheter is indicated for short-term or long-term peripheral access to the central venous system for blood sampling, intravenous administration of fluids or medications, central venous pressure monitoring and power injection of contrast media.
Target population(s)	<p>Pro-PICC®, Jet-PICC, and PFM-PICC Pro-PICC®/Jet-PICC/PFM-PICC Power Injectable Peripherally Inserted Central Catheters are intended for use in adult and pediatric patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician.</p>

	<p>Pro-PICC® Valved</p> <p>Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheters are intended for use in adult patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician. The device is not intended for use in pediatric patients.</p>
Contraindications and/or limitations	<ul style="list-style-type: none"> • The presence of device related local infection, bacteremia or septicemia is known or suspected. • The patient's body size is insufficient to accommodate the size of the implanted device. • The patient is known or is suspected to be allergic to materials contained in the device. • There has been past irradiation of prospective insertion site. • There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site. • There are local tissue factors that may prevent proper device stabilization and/or access.

3. Device description




Figure 1: Representative Image of Pro-PICC® devices

Description of device	<p>Pro-PICC®</p> <p>The Pro-PICC® Power Injectable Peripherally Inserted Central Catheter product family is available in a variety of lumen configurations and various sizes. The catheter lumen terminates at a molded hub. A proximal lumen (extension) extends from hub and terminates with a female luer-lock connector. Each extension is marked with the lumen gauge size and has a pinch clamp to control fluid flow and an ID Tag marked with the maximum power injection rate. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter. The outer diameter of the lumen increases gradually as it nears the hub. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter.</p> <p>Pro-PICC® Valved</p> <p>The Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheter family is available in a variety of lumen configurations and various sizes. The catheter lumen terminates at a molded hub. A proximal lumen</p>
-----------------------	--

	<p>(extension) extends from hub and terminates with a female luer-locking valve which controls the flow of fluids to provide clamp free infusion therapy. Positive pressure into the catheter (gravity, pump, syringe) will open the valve. When negative pressure (aspiration) is applied, the valve opens allowing the withdrawal of blood into a syringe. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter.</p> <p>Jet-PICC The Jet-PICC Power Injectable Peripherally Inserted Central Catheter product family is available in a variety of lumen configurations and various sizes. The catheter lumen terminates at a molded hub. A proximal lumen (extension) extends from hub and terminates with a female luer-lock connector. Each extension is marked with the lumen gauge size and has a pinch clamp to control fluid flow and an ID Tag marked with the maximum power injection rate. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter. The outer diameter of the lumen increases gradually as it nears the hub. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter.</p> <p>PFM-PICC The PFM-PICC Power Injectable Peripherally Inserted Central Catheter product family is available in a variety of lumen configurations and various sizes. The catheter lumen terminates at a molded hub. A proximal lumen (extension) extends from hub and terminates with a female luer-lock connector. Each extension is marked with the lumen gauge size and has a pinch clamp to control fluid flow and an ID Tag marked with the maximum power injection rate. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter. The outer diameter of the lumen increases gradually as it nears the hub. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter.</p>										
Materials / substances in contact with patient tissue	<p>The percentage ranges in the table below are based on the weight of the 3F Single Lumen (2.56g) and 6F Triple Lumen (7.45g) Pro-PICC® Power Injectable PICCs.</p> <p style="text-align: center;">Pro-PICC® Power Injectable PICCs (Non-Valved)</p> <table border="1" data-bbox="456 1482 1425 1787"> <thead> <tr> <th>Material</th><th>% Weight (w/w)</th></tr> </thead> <tbody> <tr> <td>Polyurethane</td><td>58.88 - 64.09</td></tr> <tr> <td>Acetal Co-polymer</td><td>16.82 - 24.41</td></tr> <tr> <td>Acrylonitrile Butadiene Styrene</td><td>8.13 - 10.57</td></tr> <tr> <td>Barium Sulfate</td><td>2.82 - 11.80</td></tr> </tbody> </table>	Material	% Weight (w/w)	Polyurethane	58.88 - 64.09	Acetal Co-polymer	16.82 - 24.41	Acrylonitrile Butadiene Styrene	8.13 - 10.57	Barium Sulfate	2.82 - 11.80
Material	% Weight (w/w)										
Polyurethane	58.88 - 64.09										
Acetal Co-polymer	16.82 - 24.41										
Acrylonitrile Butadiene Styrene	8.13 - 10.57										
Barium Sulfate	2.82 - 11.80										

	<p>The percentage ranges in the table below are based on the weight of the 4F Single Lumen (3.17g) and 6F Triple Lumen (7.26g) Pro-PICC® Valved Power Injectable PICCs.</p> <table><tr><th colspan="2">Pro-PICC® Power Injectable PICCs (Valved)</th></tr><tr><th>Material</th><th>% Weight (w/w)</th></tr><tr><td>Polyurethane</td><td>86.58 - 91.51</td></tr><tr><td>Barium Sulfate</td><td>5.78 - 6.83</td></tr><tr><td>Silicone</td><td>2.64 - 6.83</td></tr></table> <p>Note: Accessories containing stainless steel may contain up to 0.4% weight of the CMR substance cobalt.</p> <p>Note: Per the instructions for use, the device is contraindicated for patients with known or suspected allergies to the above materials.</p>		Pro-PICC® Power Injectable PICCs (Valved)		Material	% Weight (w/w)	Polyurethane	86.58 - 91.51	Barium Sulfate	5.78 - 6.83	Silicone	2.64 - 6.83
Pro-PICC® Power Injectable PICCs (Valved)												
Material	% Weight (w/w)											
Polyurethane	86.58 - 91.51											
Barium Sulfate	5.78 - 6.83											
Silicone	2.64 - 6.83											
Information on medicinal substances in the device	N/A											
How the device achieves its intended mode of action	<p>The subject devices utilize a Seldinger or Modified Seldinger technique to gain access. The main difference is one technique utilizes an Introduction Sheath and one does not. The Seldinger techniques for venous access are well-known surgical techniques used for inserting PICC devices. The instructions for use of each catheter are detailed in the IFUs. Catheters are to be inserted, manipulated and removed by a qualified, licensed physician or other qualified health care professional utilizing strict aseptic technique.</p> <p>Once in place, fluids are delivered or blood is withdrawn via the PICC catheter most commonly with a disposable tubing set or syringe. Catheter care includes use of a locking solution to maintain catheter patency. Catheter removal is normally done by gently pulling on the catheter, but removal may require that a surgical procedure be performed by a physician familiar with the appropriate techniques in some circumstances.</p>											
Sterilization Information	Contents sterile and non-pyrogenic in unopened, undamaged package. Sterilized by Ethylene Oxide.											
Previous generations / variants	Name of previous generation	Differences from current device										
	N/A	N/A										
Accessories intended for use in combination with the device	Name of Accessory											
	Description of Accessory											
	Part Number	Description										
	30415-018-070	0.47mm x 70cm (.018) Coated Guidewire Floppy Straight Tip										
	10129	0.76mm (0.030") I.D. Adaptor w/Sideport										
	30205-210	0.9mm OD x 0.5mm ID x 70mm (21GA) Needle W/Echo Tip										
	30824	Securement Device										
30479	Scalpel											

	30198-075	Stylet
	10700-10-035	1.1mm ID x 10cm (3.5F) Peelable Introducer
	10700-10-045	1.5mm ID x 10cm (4.5F) Peelable Introducer
	10700-10-055	1.8mm ID x 10cm (5.5F) Peelable Introducer
	10590-10-065	2.0mm ID x 10cm (6.5F) Peelable Introducer
	3035	Syringe
	3418	Tape Measure
	30823	Needleless Connector
	30415-018-13065	0.47mm x 130cm (.018) Coated Guidewire Floppy Straight Tip
	30330-018	0.47mm x 45cm (.018) Guidewire Floppy Straight Tip
	30318-021-007	0.9mm OD x 0.5mm ID x 70mm (21GA) Safety Needle W/Echo Tip
	10700-07-035	1.1mm ID x 7cm (3.5F) Peelable Introducer
	10700-07-045	1.5mm ID x 7cm (4.5F) Peelable Introducer
	10700-07-055	1.8mm ID x 7cm (5.5F) Peelable Introducer
	10590-07-065	2.0mm ID x 7cm (6.5F) Peelable Introducer

4. Risks and warnings

Residual risks and undesirable effects	As per product IFUs, All surgical procedures carry risk. Medcomp has implemented risk management processes to proactively find and mitigate these risks as far as possible without adversely affecting the benefit-risk profile of the device. After mitigation, residual risks and the possibility of adverse events from use of this product remain. Medcomp has determined that all residual risks are acceptable.	
	Residual Harm Type	Possible Adverse Events Associated with Harm
	Allergic Reaction	Allergic Reaction Intolerance Reaction to Implanted Device
	Bleeding	Bleeding Hematoma
	Cardiac Event	Cardiac Arrhythmia Cardiac Tamponade Myocardial Erosion
	Embolism	Air Embolism Thromboembolism Catheter Embolism Catheter Occlusion
	Infection	Catheter Related Sepsis Endocarditis Exit Site Infection Phlebitis
	Perforation	Perforation of Vessels or Viscus Vessel Erosion Laceration of the Vessels or Viscus
	Stenosis	Venous Stenosis

	Tissue Injury	Brachial Plexus Injury Exit Site Necrosis Soft Tissue Injury		
	Thrombosis	Venous Thrombosis Ventricular Thrombosis Fibrin Sheath Formation		
	Miscellaneous complications	Catheter Erosion Through Skin Spontaneous Catheter Tip Malposition or Retraction Risks Normally Associated with Local or General Anesthesia, Surgery and Post-Operative Recovery		
	Patient Residual Harm Category	Quantification of Residual Risks		
		PMS Complaints (01 January 2019 – 30 September 2024)	PMCF Events	
		Units Sold: 670,138	Units Studied: 2,030	
		% of Devices	% of Devices	
		Allergic Reaction	Not Reported	Not Reported
		Bleeding	0.00015%	Not Reported
		Cardiac Event	Not Reported	Not Reported
		Embolism	Not Reported	Not Reported
		Infection	Not Reported	1.77%
		Perforation	Not Reported	Not Reported
Stenosis		Not Reported	Not Reported	
Tissue Injury		0.00015%	Not Reported	
Thrombosis	Not Reported	Not Reported		
Warnings and precautions	All warnings have been reviewed against the risk analysis, PMS, and usability testing to validate consistency between the sources of information.			
	<ul style="list-style-type: none">Do not insert catheter in thrombosed vesselsDo not advance the guidewire or catheter if unusual resistance is encountered.Do not insert or withdraw the guidewire forcibly from any component. If the guidewire becomes damaged, guidewire and any associated componentry must be removed together.Do not resterilize the catheter or accessories by any method.Contents sterile and non-pyrogenic in unopened, undamaged package. STERILIZED BY ETHYLENE OXIDEDo not re-use catheter or accessories as there may be a failure to adequately clean and decontaminate the device which may lead to contamination, catheter degradation, device fatigue, or endotoxin reaction.Do not use catheter or accessories if package is opened or damaged.			

- Do not use catheter or accessories if any sign of product damage is visible or the use-by date has passed.
- Do not use sharp instruments near the extension lines or catheter lumen.
- Do not use scissors to remove dressing.

Precautions listed in the IFUs are as follows:

- The fluid level of the catheter will drop (allowing air entry) if the catheter connector is held above the level of the patient's heart and opened to air. To help prevent a drop in the fluid volume (allowing air entry) while changing injection caps, hold the connector below the level of the patient's heart before removing the injection cap.*
- Refer to standards of practice and institutional policies for compatible infusion agents for central venous access.
- Follow all contraindications, warnings, precautions, and instructions for all infusates as specified by their manufacturer.
- Syringes smaller than ten (10) ml will generate excessive pressure and may damage the catheter. Ten (10) ml or larger syringes are recommended.
- Hydrate guidewire prior to use.
- Always flush catheter prior to removing stylet.**
- Always flush catheter prior to insertion.*
- Catheter will be damaged if clamps other than what is provided with this kit are used.**
- Catheter will be damaged if clamps are used.*
- Clamping of the tubing repeatedly in the same location may weaken tubing. Avoid clamping near the luer(s) and hub of the catheter.**
- Examine catheter lumen and extension(s) before and after each infusion for damage.
- To prevent accidents, assure the security of all caps and connections prior to and between uses.
- Use only Luer Lock (threaded) Connectors with this catheter.
- In the rare event that a hub or connector separates from any component during insertion or use, take all necessary steps and precautions to prevent blood loss or air embolism.
- Repeated overtightening of bloodlines, syringes, and caps will reduce connector life and could lead to potential connector failure.
- Confirm catheter tip position prior to use. Monitor tip placement routinely per institution policy.
- This is not a right atrium catheter. Avoid positioning the catheter tip in the right atrium. Placement or migration of the catheter tip into the right atrium may cause cardiac arrhythmia, myocardia erosion, or cardiac tamponade.
- The valve is not a barrier to infection. Strict aseptic technique must be utilized during all actuations and cap changes. A sterile

	<p>end cap should be applied to the hub of the catheter to prevent contamination when not in use.*</p> <ul style="list-style-type: none"> • Discard biohazard according to facility protocol. • The CMR substance Cobalt is a naturally occurring component of stainless steel. Based on biocompatibility evaluation it was determined that the main hazards of stainless steels are related to the processing of the material, especially welding, thus not applicable to the intended use of the device. Stainless steels used in these devices are unlikely to reach exposure levels that will elicit carcinogenicity, mutagenicity or reproductive toxicity. <p><i>*Precaution only appears in Pro-PICC® Valved IFU (40798BSI)</i> <i>**Precaution only appears in Pro-PICC® IFUs (40795BSI, 40795JBSI, 40795PBSI)</i></p>
Other relevant aspects of safety (ex. field safety corrective actions, etc.)	<p>For a period of 01 January 2019 to 30 September 2024 there were 186 complaints for 670,138 units sold, giving an overall complaint rate of 0.028%. There were no death-related events. No events resulted in recalls during the review period.</p>

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)

Summary of clinical data related to the subject device			
Clinical Literature	PMCF Data	Total	User Survey Responses
467 (& 3,580 Mixed Cohort Cases)	2,030	2,497 (& 3,580 Mixed Cohort Cases)	37
<p>Clinical performance was measured using parameters including but not limited to dwell time, catheter insertion outcomes, and adverse event rates. Critical clinical parameters extracted from these studies met standards set forth in the guidelines for the State of the Art. There were no unforeseen adverse events or other high occurrences of adverse events detected in any of the clinical activities.</p> <p>Survivability of a given implant is a multi-factorial event that depends on numerous factors, including: the limits of the implant, surgical technique, difficulty level of the surgical procedure, patient health, patient activity level, patient medical history, and other factors. In the case of the Pro-PICC® Power Injectable Peripherally Inserted Central Catheter, 93 catheters had a 55.07 day [95%CI: 43.98-66.18 days] duration of use that has been found in clinical use reported to date. Based on this information, the Pro-PICC® Power Injectable Peripherally Inserted Central Catheter has a 12 month lifetime; however, the decision to remove and/or replace the catheter should be based on clinical performance and need, and not any predetermined point in time.</p>			
Summary of clinical data related to the equivalent device (if applicable)			
<p>Clinical evidence from published literature and PMCF activities has been generated specific to known and unknown variants of the subject device. The equivalency rationale in the</p>			

updated clinical evaluation report will demonstrate that the clinical evidence available for these variants is representative of the range of device variants in the device family.

There are no clinical or biological differences between variants within the subject device family, and the potential impact of the technical differences will be rationalized in the updated clinical evaluation report.

Summary of clinical data from pre-market investigations (if applicable)

No pre-market clinical investigations were used for the device's clinical evaluation.

Summary of clinical data from other sources:

Source: Summary of Published Literature

Clinical evidence literature searches have found sixteen published literature articles representing 467 Pro-PICC® device family specific cases, and an additional 3,580 mixed cohort cases inclusive of the Pro-PICC® device family. The articles included randomized prospective studies (Paquet et al., Pittiruti et al., Yong et al.), prospective studies (Cotogni et al., Derudas et al., Zerla et al.), retrospective studies (Annetta et al., Annetta et al., Jeon et al., Kim et al., Wortley et al., Yang et al., Yeon et al., Yu and Hong), a hybrid prospective/retrospective study (Biasucci et al.), and a conference proceeding (Casas et al.).

Bibliography:

- Annetta MG, Marche B, Dolcetti L, et al. Ultrasound-guided cannulation of the superficial femoral vein for central venous access. *J Vasc Access*. 2021
- Annetta MG, Ostroff M, Marche B, et al. Chest-to-arm tunneling: A novel technique for medium/long term venous access devices. *J Vasc Access*. 2021
- Biasucci, D. G., Pittiruti, M., Taddei, A., Picconi, E., Pizza, A., Celentano, D., . . . Conti, G. (2018). Targeting zero catheter-related bloodstream infections in pediatric intensive care unit: a retrospective matched case-control study. *The journal of vascular access*, 19(2), 119-124.
- Cotogni, P., Barbero, C., Garrino, C., Degiorgis, C., Mussa, B., De Francesco, A., & Pittiruti, M. (2015). Peripherally inserted central catheters in non-hospitalized cancer patients: 5-year results of a prospective study. *Support Care Cancer*, 23(2), 403-409. doi:10.1007/s00520-014-2387-9
- Cotogni P, Mussa B, Degiorgis C, De Francesco A, Pittiruti M. Comparative Complication Rates of 854 Central Venous Access Devices for Home Parenteral Nutrition in Cancer Patients: A Prospective Study of Over 169,000 Catheter-Days. *Journal of Parenteral and Enteral Nutrition*. 2021;45(4):768-76.
- Derudas, D., Massidda, S., Simula, M. P., Dessì, D., Usai, S. V., Longhitano, G., Daniela Ibba, Loredana Aracu, Monica Atzori & La Nasa, G. (2023). Peripherally inserted central catheter insertion and management in Hodgkin and non-Hodgkin lymphomas: a 13-year monocentric experience. *Frontiers in Hematology*, 2, 1171991.
- Jeon, E.-Y., Cho, Y. K., Yoon, D. Y., & Hwang, J. H. (2016). Which arm and vein are more appropriate for single-step, non-fluoroscopic, peripherally inserted central catheter insertion? *The journal of vascular access*, 17(3), 249-255.
- Kim, H., Cho, S. B., Park, S. E., Jo, S. H., Lim, S. G., Jeong, Y., JH Won, WJ Yang, HC Choi, JH Ahn & Nam, I. C. (2024). A New Equation to Estimate Peripherally Inserted Central Catheter Length. *Medicina*, 60(3), 417.

- Paquet, F., Boucher, L. M., Valenti, D., & Lindsay, R. (2017). Impact of arm selection on the incidence of PICC complications: results of a randomized controlled trial. *J Vasc Access*, 18(5), 408-414. doi:10.5301/jva.5000738
- Pittiruti, M., Emoli, A., Porta, P., Marche, B., DeAngelis, R., & Scoppettuolo, G. (2014). A prospective, randomized comparison of three different types of valved and non-valved peripherally inserted central catheters. *J Vasc Access*, 15(6), 519-523. doi:10.5301/jva.5000280
- Sze Yong T, Vijayanathan AA, Chung E, Ng WL, Yaakup NA, Sulaiman N. Comparing catheter related bloodstream infection rate between cuffed tunnelled and non-cuffed tunnelled peripherally inserted central catheter. *J Vasc Access*. 2022;23(2):225-31.
- Wortley, V., & Almerol, L. A. (2020). Misplacement of piccs following power-injected CT contrast media. *British Journal of Nursing*, 29(19), S4-S10.
- Yang WJ, Kang D, Shin JH, et al. Comparison of different techniques for the management of venous steno-occlusive lesions during placement of peripherally inserted central catheter. *Sci Rep*. 2021;11(1)
- Yeon, J. W., Cho, Y. K., Kim, H. M., Song, M. G., Song, S.-Y., Cho, S. B., & Lee, S. Y. (2018). Interventional management of central vein occlusion in patients with peripherally inserted central catheter placement. *Journal of Vascular Surgery: Venous and Lymphatic Disorders*, 6(5), 566-574.
- Yu, B., & Hong, J. (2022). Safety and Efficacy of Peripherally Inserted Central Catheter Placement by Surgical Intensivist–Led Vascular Access Team. *Vascular Specialist International*, 38(4).
- Zerla, P. A., Canelli, A., Cerne, L., Caravella, G., Gilardini, A., De Luca, G., . . . Venezia, R. (2017). Evaluating safety, efficacy, and cost-effectiveness of PICC securement by subcutaneously anchored stabilization device. *The journal of vascular access*, 18(3), 238-242.

Source: PMCF_Infusion_201

The CVAD Registry was acquired from CVAD Resources, LLC on 23 August 2020. All data received was de-identified, but otherwise represented exactly what was entered by clinicians on a consecutive basis. Medcomp received only data pertaining to devices with the manufacturer listed as “Medcomp” and all case information was sourced from two US hospitals. Hospital ID 121 is described as a “Vascular Access team in a Not-for-Profit Community Based Hospital”, and Hospital ID 123 is described as a “PICC (peripherally inserted central catheter) team in an Academic Medical Center”. Insertion of device dates range from 06 August 2012 through 21 April 2015. Removal of device dates from 09 August 2012 through 07 May 2015.

1,826 Pro-PICC® cases inclusive of several variant devices across French size (4F, 5F, and 6F) and lumen configuration (single, double, triple) were collected. The following outcome measures were confirmed to be within State of the Art safety and performance outcome measures from published literature for Medcomp Pro-PICC® devices:

- Dwell Time – 13.5 Days (95%CI: 11.8 – 15.2)
- Procedural Outcomes – 98.6% (95%CI: 98.1% - 99.1%)
- Catheter Related Blood Stream Infection – 2.4 Confirmed Events per 1,000 Catheter Days

The variants included in the dataset are displayed below.

Variant	n	French sizes
Single Lumen Pro-PICC	30	4F, 5F, 6F
Double Lumen Pro-PICC	1647	5F, 6F
Triple Lumen Pro-PICC	129	6F
Pro-PICC Unknown	20	5F

Source: PMCF_Infusion_211

The Infusion Product Line Data Collection Survey aimed to assess safety and performance outcome information for all variants of Medcomp Infusion Ports, PICCs, Midlines, and CVCs. 70 survey responses were collected from 17 countries representing 471 device cases.

204 Pro-PICC® cases inclusive of several variant devices across French size (3F, 4F, 5F, and 6F), Valve (with and without), and lumen configuration (single, double, triple) were collected. The following outcome measures were confirmed to be within State of the Art safety and performance outcome measures from published literature for Medcomp Pro-PICC® devices:

- Dwell Time – 55.07 Days (95%CI: 43.98 – 66.18)
- Procedural Outcomes – 95.10% (95%CI: 94.4% - 95.8%)
- Phlebitis – No Events Reported
- Infiltration/Extravasation – No Events Reported
- Catheter Associated Venous Thrombus – No Events Reported
- Catheter Related Blood Stream Infection – 0.39 per 1,000 catheter days (95%CI: 0 – 0.93)
- Power Injection Related Complications – 0.43% (95%CI: 0% - 1.3%)

The variants included in the dataset are displayed below.

Variant	n	French Size(s)	Length(s)
Single Lumen Pro-PICC	105	3F, 4F, 5F	55cm, 60cm
Single Lumen Valved Pro-PICC	5	5F	60cm
Dual Lumen Pro-PICC	68	4F, 5F, 6F	55cm, 60cm
Triple Lumen Pro-PICC	26	6F	60cm

Source: Duration of Use Customer Survey

An email questionnaire was distributed globally to users of Medcomp PICCs and CVCs from 10 October 2019 to 16 October 2019. The questionnaire asked respondents to identify, from their own experience, the number of products used yearly, the average dwell time, and the longest dwell time for each applicable device family.

Across the five device families, a total of 69 responses were collected from 14 countries. Means and ranges of responses for each device family were compiled on 16 October 2019.

24 responses were acquired relating to the Pro-PICC® device family. Across an estimated 8,761 products used yearly, the mean average dwell time was 116 days (Range: 14 – 365 days), and the mean longest dwell time was 360 days (Range: 60 – 2,555 days).

Source: PMCF_Medcomp_211

The Medcomp User Survey acquired responses from healthcare personnel familiar with any number of Medcomp's product offerings.

13 respondents responded that they or their facility have used Medcomp PICCs, with 13 of those respondents using the Pro-PICC® device. There were no differences in mean user sentiments within PICCs across State of the Art Performance and Safety Outcome Measures or between device types relating to safety or performance.

The following data points were collected from users of Medcomp PICCs (n=13):

- (Mean Likert Scale Response) Catheters function as intended – 4.7 / 5
- (Mean Likert Scale Response) Packaging allows for aseptic presentation – 4.9 / 5
- (Mean Likert Scale Response) Benefit outweighs the risk – 4.6 / 5
- Dwell Time (n=11) – 58.1 days (**95%CI:** 15.5 – 100.8)

The following data points were collected from users of Medcomp Pro-PICC® (n=11):

- (Mean Likert Scale Response) Catheters function as intended – 4.7 / 5
- (Mean Likert Scale Response) Packaging allows for aseptic presentation – 4.9 / 5
- (Mean Likert Scale Response) Benefit outweighs the risk – 4.7 / 5
- Dwell Time (n=8) – 66 days (**95%CI:** 3.7 – 128.3)

The following complications were reported for Pro-PICC® devices:

- Placement Issues (No Comments on Frequency)
- DVT (Deep Vein Thrombosis) (No Comments on Frequency)
- Infection (No Comments on Frequency)
- Thrombosis (No Comments on Frequency)
- Occlusion (No Comments on Frequency)

Overall summary of clinical safety and performance

Upon review of the data across all sources, it is possible to conclude that the benefits of the subject device, which is facilitating blood withdraws for laboratory testing, delivery of fluids and medications for treatments including chemotherapy, and power injection of contrast media for CT examinations in patients in whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks deemed necessary based on the direction of a qualified, licensed physician, outweigh the overall and individual risks when the device is used as intended by the manufacturer. It is the manufacturer's and clinical expert evaluator's opinion that activities both complete and ongoing are sufficient to support the safety, efficacy, and acceptable benefit/risk profile of the subject devices.

Outcome	Benefit/Risk Acceptability Criteria	Desired Trend	Clinical Literature (Subject Device)	PMCF Data (Subject Device)
Performance				
Dwell Time	Greater than 6.27 days	↑	15 – 575 Days (Summary)	55.07 Days (PMCF_Infusion_211)

			of Published Literature)	13.5 Days (PMCF_Infusion_201) 116 Days (Duration of Use Customer Survey) 66 Days (PMCF_Medcomp_211) Likert Scale Response 4.7 / 5 (Section 6.5.8)**
Procedural Outcomes	Greater than 43% (Bedside) / 90% (Interventional Radiology)	↑	46.3% - 53.7% at Bedside (Summary of Published Literature) 99.7% - 100% (Summary of Published Literature)	95.10% (PMCF_Infusion_211) 98.6% (PMCF_Infusion_201) Likert Scale Response 4.7 / 5 (PMCF_Medcomp_211)**
Safety				
Phlebitis	Less than 2.4% catheters with reported incidents of phlebitis	↓	ND*	None Reported (PMCF_Infusion_211) Likert Scale Response 4.6 / 5 (PMCF_Medcomp_211)**
Infiltration/Extravasation	Less than 7% catheters with reported incidents of infiltration or extravasation	↓	0.6% - 7% (Summary of Published Literature)	None Reported (PMCF_Infusion_211) Likert Scale Response 4.7 / 5 (PMCF_Medcomp_211)**
Catheter Associated Venous Thrombus (CAVT)	Less than 5.4 incidents of CAVT per 1,000 catheter days	↓	0 – 1.0 per 1,000 catheter days (Summary of Published Literature)	None Reported (PMCF_Infusion_211) Likert Scale Response 4.5 / 5 (PMCF_Medcomp_211)**
Central Line Associated Blood Stream Infection (CLABSI) / Catheter Related Blood Stream Infection (CRBSI)	Less than 5.7 incidents of CLABSI/CRBSI per 1,000 catheter days	↓	0 – 2.73 per 1,000 catheter days (Summary of Published Literature)	0.39 per 1,000 catheter days (PMCF_Infusion_211) 2.4 per 1,000 catheter days (PMCF_Infusion_201)

				Likert Scale Response 4.7 / 5 (PMCF_Medcomp_211)**
Power Injection Related Complications	Less than 1.8% catheters with reported incidents of rupture due to contrast injection Less than 15.4% catheters with reported incidents of displacement due to contrast injection	↓	0.6% - 0.7% (Summary of Published Literature)	0.43% (PMCF_Infusion_211) Likert Scale Response 4.7 / 5 (PMCF_Medcomp_211)**

*ND indicates no data on the clinical data parameter

**PMCF_Medcomp_211 asked respondents, if they agreed on a scale of 1 -5, that their experience in relation to each outcome was the same or better than the benefit/risk acceptability criteria.

On-going or planned Post-Market Clinical Follow-up (PMCF)

Activity	Description	Reference	Timeline
Multicenter Patient-Level Case Series	Collect additional clinical data on the device	PMCF_PICC_231	Q4 2025
State of the Art Literature Search	Identify risks and trends with use of similar devices	SAP-Infusion	Q2 2025
Clinical Evidence Literature Search	Identify risks and trends with use of the device	LRP-Infusion	Q2 2025
Global Trial Database Search	Identify ongoing clinical trials involving Medcomp® catheters	N/A	Q3 2025
Truveta Data Queries and Retrospective Analysis	Collect additional clinical data on the device and comparators	TBD	Q4 2025

No emerging risks, complications or unexpected device failures have been detected from PMCF activities.

6. Possible therapeutic alternatives

The Infusion Nurses Society (INS) Standards 2021 clinical practice guidelines have been used to support the below recommendations for treatments.

Therapy	Benefits	Disadvantages	Key Risks
Central Venous Catheters (CVCs)	<ul style="list-style-type: none"> Easy access once in place Minimizes repeated venipuncture Increased patient mobility during infusion 	<ul style="list-style-type: none"> Requires surgical procedure for placement Risks associated with surgery: general anesthesia, etc. Requires maintenance 	<ul style="list-style-type: none"> Catheter infection <ul style="list-style-type: none"> Occlusion Malfunction of the CVC Vascular thrombosis

Therapy	Benefits	Disadvantages	Key Risks
	<ul style="list-style-type: none"> Easier for outpatient treatment 	<ul style="list-style-type: none"> High risk of infection or thrombotic event 	
Implantable Ports	<ul style="list-style-type: none"> Decreases puncture wounds/vein damage compared to traditional injection Easier to visualize, palpate, and therefore safer form of IV access Reduces chance for corrosive medications to make skin contact <ul style="list-style-type: none"> Only one venipuncture for both treatment and lab draws, as opposed to two for traditional IV Longer dwelling time compared to IV <ul style="list-style-type: none"> Can be permanent, if needed 	<ul style="list-style-type: none"> Requires surgical procedure, but IV does not Risks associated with surgery: general anesthesia, etc. Requires regular flushing 	<ul style="list-style-type: none"> Drug extravasations <ul style="list-style-type: none"> Infection Thromboembolism Tissue necrosis of overlying skin / port dehiscence
Midline Catheters	<ul style="list-style-type: none"> Patient comfort – fewer restarts than IVs Longer dwell time than IVs Lower risk of infection compared to IVs <ul style="list-style-type: none"> No X-ray required before use Decreased chance of extravasation of infusate 	<ul style="list-style-type: none"> Data on clear disadvantages compared to other modalities is not available Not suitable for continuous injections of most vesicants or irritants 	<ul style="list-style-type: none"> Insertion-related phlebitis
Peripherally Inserted Central	<ul style="list-style-type: none"> Decreased risk of catheter occlusion compared to CVC 	<ul style="list-style-type: none"> Increased risk of deep vein thrombosis compared to CVC 	<ul style="list-style-type: none"> Deep vein thrombosis (DVT)

Therapy	Benefits	Disadvantages	Key Risks
Catheters (PICCs)	<ul style="list-style-type: none"> Fewer venous punctures compared to traditional PIV 	<ul style="list-style-type: none"> Pain/Discomfort over time Need for adaptation in daily life 	<ul style="list-style-type: none"> Pulmonary embolism Venous thromboembolism (VTE) Post thrombotic syndrome
Peripheral Intravenous Catheters (PIVs)	<ul style="list-style-type: none"> Does not require surgical procedure 	<ul style="list-style-type: none"> Higher hemolysis rates compared to venipuncture Infection Hematoma/thrombosis Cannot be used for therapies with blistering agents Four days maximum use 	<ul style="list-style-type: none"> Infection Phlebitis

7. Suggested profile and training for users

The catheter should be inserted, manipulated, and removed by a qualified, licensed physician or other qualified health care professional under the direction of a physician.

8. Reference to any harmonized standards and Common Specifications (CS) applied

Harmonized Standard or CS	Revision	Title or Description	Level of Compliance
EN 556-1	2001	Sterilization of medical devices. Requirements for medical devices to be designated "STERILE". Requirements for terminally sterilized medical devices	Full
EN ISO 10555-1	2013+A1:2017	Intravascular catheters. Sterile and single-use catheters. General requirements	Full
EN ISO 10555-3	2013	Intravascular catheters. Sterile and single-use catheters. Central venous catheters	Full
EN ISO 10993-1	2020	Biological evaluation of medical devices — Part 1: Evaluation and testing within a risk management process	Full
EN ISO 10993-7	2008+ A1:2022	Biological evaluation of medical devices — Part 7: Ethylene oxide sterilization residuals — Amendment 1: Applicability of allowable limits for neonates and infants	Full
EN ISO 10993-18	2020	Biological evaluation of medical devices — Part 18: Chemical characterization of medical device materials within a risk management process	Full
EN ISO 11070	2014+A1:2018	Sterile single-use intravascular introducers, dilators and guidewires	Full
EN ISO 11135	2014 + A1: 2019	Sterilization of health-care products. Ethylene oxide. Requirements for the development, validation and routine	Full

Harmonized Standard or CS	Revision	Title or Description	Level of Compliance
		control of a sterilization process for medical devices	
EN ISO 11138-1	2017	Sterilization of health care products — Biological indicators Part 1: General requirements	Full
EN ISO 11138-2	2017	Sterilization of health care products— Biological indicators—Part 2: Biological indicators for ethylene oxide sterilization processes	Full
EN ISO 11138-7	2019	Sterilization of health care products. Biological indicators - Guidance for the selection, use and interpretation of results	Full
EN ISO 11140-1	2014	Sterilization of health care products — Chemical indicators Part 1: General requirements	Full
EN ISO 11607-1 Excludes Section 7	2020	Packaging for terminally sterilized medical devices. Requirements for materials, sterile barrier systems and packaging systems	Partial; (Transition Plan)
EN ISO 11607-2	2020	Packaging for terminally sterilized medical devices. Validation requirements for forming, sealing and assembly processes	Full
EN ISO 11737-1	2018 + A1: 2021	Sterilization of health care products. Microbiological methods. Determination of a population of microorganisms on products	Full
EN ISO 13485	2016 + A11: 2021	Medical Devices – Quality Management system – Requirements for Regulatory Purposes	Full
EN ISO 14155	2020	Clinical investigation of medical devices for human subjects — Good clinical practice	Full
EN ISO 14644-1	2015	Cleanrooms and associated controlled environments — Part 1: Classification of air cleanliness by particle concentration	Full
EN ISO 14644-2	2015	Cleanrooms and associated controlled environments — Part 2: Monitoring to provide evidence of cleanroom performance related to air cleanliness by particle concentration	Full
EN ISO 14971	2019+A11:2021	Medical devices. Application of risk management to medical devices	Full
EN ISO 15223-1	2021	Medical devices — Symbols to be used with medical device labels, labelling and information to be supplied — Part 1: General requirements	Full
EN ISO/IEC 17025	2017	General requirements for the competence of testing and calibration laboratories	Full
PD CEN ISO/TR 20416	2020	Medical devices — post-market surveillance for manufacturers	Full
EN ISO 20417	2021	Medical devices - Information to be supplied by the manufacturer.	Full

Harmonized Standard or CS	Revision	Title or Description	Level of Compliance
EN 62366-1	2015 + A1: 2020	Medical devices — Part 1: Application of usability engineering to medical devices	Full
ISO 7000	2019	Graphical symbols for use on equipment. Registered symbols	Partial
ISO 594-1	1986	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment — Part 1: General requirements	Full
ISO 594-2	1998	Conical fittings with a 6 % (Luer) taper for syringes, needles and certain other medical equipment — Part 2: Lock Fittings	Full
MEDDEV 2.7.1	Rev 4	Clinical Evaluation: A Guide for Manufacturers and Notified Bodies Under Directives 93/42/EEC and 90/385/EEC	Full
MEDDEV 2.12/2	Rev. 2	GUIDELINES ON MEDICAL DEVICES POST MARKET CLINICAL FOLLOW-UP STUDIES A GUIDE FOR MANUFACTURERS AND NOTIFIED BODIES	Full
MDCG 2020-6	2020	Clinical evidence needed for medical devices previously CE marked under Directives 93/42/EEC or 90/385/EEC	Full
MDCG 2020-7	2020	Post-market clinical follow-up (PMCF) Plan Template A guide for manufacturers and notified bodies	Full
MDCG 2020-8	2020	Post-market clinical follow-up (PMCF) Evaluation Report Template A guide for manufacturers and notified bodies	Full
MDCG 2019-9	2022	Summary of safety and clinical performance	Full
MDCG 2018-1	Rev. 4	Guidance on BASIC UDI-DI and changes to UDI-DI	Full
ASTM D 4169-16	2022	Standard Practices for Performance Testing of Shipping Containers and Systems.	Full
ASTM F2096-11	2019	Standard Test Method for Detecting Gross Leaks in Packaging by Internal Pressurization (Bubble Test)	Full
ASTM F2503-20	2020	Standard Practice for Marking Medical Devices and Other Items for Safety in the Magnetic Resonance Environment	Full
ASTM F640-20	2020	Standard Test Methods for determining Radiopacity for Medical Use	Full
ASTM D4332-14	2014	Standard Practice for Conditioning Containers, Packages, or Packaging Components for Testing	Full

PATIENTS

SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Revision: SSCP-013 Rev. 5

Date: 24OCT2024

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

IMPORTANT INFORMATION

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your healthcare professional in case you have questions about your medical condition or about the use of the device in your situation.

This SSCP is not intended to replace an Implant Card or the Instructions for Use to provide information on the safe use of the device.

1. Device identification and general information

Device trade name(s)	Pro-PICC® Power Injectable Peripherally Inserted Central Catheter
Manufacturer name and address	Medical Components, Inc. 1499 Delp Drive Harleysville, PA 19438, United States of America (USA)
Basic UDI-DI	00884908286NP
Date first CE certificate was issued for this device	Pro-PICC® - October 2007 Pro-PICC® Valved - May 2013 Jet-PICC - July 2009 PFM-PICC - December 2015

The devices in scope of this document are all peripherally inserted central catheter (PICC) sets. The catheter part numbers are organized into variant categories. These devices are distributed as procedure trays. Procedure trays come in different configurations.

Variant Devices:

Variant Description	Part Number(s)
3F x 55cm Single Lumen Pro-PICC®	10467-855-800 10467-855-801
4F x 55cm Single Lumen Pro-PICC®	10602-855-800 10602-855-801
4F x 55cm Single Lumen Valved Pro-PICC®	10643-855-801
5F x 55cm Double Lumen Pro-PICC®	10561-855-800 10561-855-801
5F x 55cm Double Lumen Valved Pro-PICC®	10645-855-801
5F x 60cm Single Lumen Pro-PICC®	10556-860-800 10556-860-801
5F x 60cm Single Lumen Valved Pro-PICC®	10644-860-801
6F x 60cm Double Lumen Pro-PICC®	10563-860-800 10563-860-801
6F x 60cm Triple Lumen Pro-PICC®	10568-860-800 10568-860-801
6F x 60cm Triple Lumen Valved Pro-PICC®	10646-860-801

Procedure Trays:

Catalog Code	Part Number	Description
JSACT5D	10561-855-800	5F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
JSACT5DL	10561-855-800	5F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
JSACT6D	10563-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
JSACT6DL	10563-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
PFMCT5DLWS	10561-855-800	5F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT5DS	10561-855-800	5F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP52024	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP52028	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MRCTP62024	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET
MRCTP62028	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MR17035201	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR17035202	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
MR17035205	10561-855-801	5F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN NURSING SET

Catalog Code	Part Number	Description
MR17036201	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR17036202	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN LONG WIRE SET
MR17036205	10563-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN NURSING SET
JSACT4S	10602-855-800	4F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
JSACT4SL	10602-855-800	4F X 55CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
JSACT5S	10556-860-800	5F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
JSACT5SL	10556-860-800	5F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
PFMCT3SS	10467-855-800	3F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
PFMCT4SLWS	10602-855-800	4F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT4SS	10602-855-800	4F X 55CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
PFMCT5SLWS	10556-860-800	5F X 60CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
PFMCT5SS	10556-860-800	5F X 60CM PFM-PICC POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP31024	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP41024	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP41028	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MRCTP51024	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET
MRCTP51028	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN FULL CUT DOWN SET WITH 70CM GUIDEWIRE
MR17033101	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR17033102	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17033105	10467-855-801	3F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
MR17034101	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR17034102	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17034105	10602-855-801	4F X 55CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
MR17035101	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET

Catalog Code	Part Number	Description
MR17035102	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN LONG WIRE SET
MR17035105	10556-860-801	5F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN NURSING SET
JSACT6T	10568-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET
JSACT6TL	10568-860-800	6F X 60CM JET-PICC CT POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN LONG WIRE SET
MRCTP63024	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN FULL CUT DOWN SET
MR17036301	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET
MR17036302	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN LONG WIRE SET
MR17036305	10568-860-801	6F X 60CM PRO-PICC® POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN NURSING SET
MR82034101	10643-855-801	4F X 55CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR82035101	10644-860-801	5F X 60CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER SINGLE LUMEN BASIC SET
MR82035201	10645-855-801	5F X 55CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER DUAL LUMEN BASIC SET
MR82036301	10646-860-801	6F X 60CM PRO-PICC® VALVED POWER INJECTABLE PERIPHERALLY INSERTED CENTRAL CATHETER TRIPLE LUMEN BASIC SET

Configurations of Procedure Trays:

Configuration Type
Pro-PICC® Basic Set
Pro-PICC® Longwire Set
Pro-PICC® Nursing Set
Pro-PICC® Full Cut Down Set
Pro-PICC® Full Cut Down Set with 70 centimeter (cm) Guidewire
Pro-PICC® Valved Basic Set
Jet-PICC CT Basic Set
Jet-PICC CT Longwire Set
PFM-PICC Full Cut Down Set
PFM-PICC Full Cut Down Set with 70cm Guidewire

2. Intended use of the device

Intended purpose	Pro-PICC®, Jet-PICC, and PFM-PICC Pro-PICC®/Jet-PICC/PFM-PICC Power Injectable Peripherally Inserted Central Catheters are intended for use in adult and pediatric patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified,
------------------	--

	<p>licensed physician. The device is intended to be used under the regular review and assessment of qualified health professionals. This catheter is for Single Use Only.</p> <p>Pro-PICC® Valved Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheters are intended for use in adult patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician. The device is intended to be used under the regular review and assessment of qualified health professionals. This catheter is for Single Use Only.</p>
Indication(s)	<p>The Pro-PICC®/Pro-PICC® Valved/PFM-PICC/Jet-PICC Power Injectable Peripherally Inserted Central Catheter is indicated for short-term or long-term peripheral access to the central venous system for blood sampling, intravenous administration of fluids or medications, central venous pressure monitoring and power injection of contrast media.</p>
Intended patient group(s)	<p>Pro-PICC®, Jet-PICC, and PFM-PICC Pro-PICC®/Jet-PICC/PFM-PICC Power Injectable Peripherally Inserted Central Catheters are intended for use in adult and pediatric patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician.</p> <p>Pro-PICC® Valved Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheters are intended for use in adult patients requiring frequent needlesticks for whom short-term or long-term peripheral access to the central venous system without requiring frequent needlesticks is deemed necessary based on the direction of a qualified, licensed physician. The device is not intended for use in pediatric patients.</p>
Contraindications	<ul style="list-style-type: none"> • The presence of device related local infection, bacteremia or septicemia is known or suspected. • The patient's body size is insufficient to accommodate the size of the implanted device. • The patient is known or is suspected to be allergic to materials contained in the device. • There has been past irradiation of prospective insertion site. • There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site. • There are local tissue factors that may prevent proper device stabilization and/or access.

3. Device description




Figure 1: Pro-PICC® Power Injectable Peripherally Inserted Central Catheter

Description of device	<p>Pro-PICC®</p> <p>The Pro-PICC® Power Injectable Peripherally Inserted Central Catheter product family is available in many of lumen shapes and various sizes. The catheter lumen ends at a molded hub. A proximal lumen (extension) extends from hub and ends with a female connector. Each extension is marked with the lumen gauge size and has a pinch clamp to control fluid flow and an identification (ID) Tag marked with the highest power injection rate. The outer diameter of the lumen gets bigger as it nears the hub. The highest recommended infusion rate varies by catheter French size and is printed on the catheter. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter.</p> <p>Pro-PICC® Valved</p> <p>The Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheter family is available in a variety of lumen shapes and various sizes. The catheter lumen ends at a molded hub. A proximal lumen (extension) extends from hub and ends with a female valve which controls the flow of fluids to provide clamp free infusion therapy. Positive pressure into the catheter (gravity, pump, syringe) will open the valve. When negative pressure (aspiration) is applied, the valve opens as allows blood to be drawn into a syringe. The highest recommended infusion rate varies by catheter French size and is printed on the catheter.</p> <p>Jet-PICC</p> <p>The Jet-PICC Power Injectable Peripherally Inserted Central Catheter product family is available in a variety of lumen shapes and various sizes. The catheter lumen ends at a molded hub. A proximal lumen (extension) extends from hub and ends with a female connector. Each extension is marked with the lumen gauge size and has a pinch clamp</p>
-----------------------	--

	<p>to control fluid flow and an ID Tag marked with the maximum power injection rate. The outer diameter of the lumen increases gradually as it nears the hub. The highest recommended infusion rate varies by catheter French size and is printed on the catheter. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter</p> <p>PFM-PICC The PFM-PICC Power Injectable Peripherally Inserted Central Catheter product family is available in a variety of lumen shapes and various sizes. The catheter lumen ends at a molded hub. A proximal lumen (extension) extends from hub and terminates with a female connector. Each extension is marked with the lumen gauge size and has a pinch clamp to control fluid flow and an ID Tag marked with the maximum power injection rate. The outer diameter of the lumen increases gradually as it nears the hub. The highest recommended infusion rate varies by catheter French size and is printed on the catheter. The lumen is marked with depth marks every centimeter and numerical markings every fifth centimeter.</p>																		
Materials / substances in contact with patient tissue	<p>The percentage ranges in the table below are based on the weight of the 3F Single Lumen (2.56g) and 6F Triple Lumen (7.45g) Pro-PICC® Power Injectable PICCs.</p> <p style="text-align: center;">Pro-PICC® Power Injectable PICCs (Non-Valved)</p> <table border="1"> <thead> <tr> <th>Material</th><th>% Weight (w/w)</th></tr> </thead> <tbody> <tr> <td>Polyurethane</td><td>58.88 - 64.09</td></tr> <tr> <td>Acetal Co-polymer</td><td>16.82 - 24.41</td></tr> <tr> <td>Acrylonitrile Butadiene Styrene</td><td>8.13 - 10.57</td></tr> <tr> <td>Barium Sulfate</td><td>2.82 - 11.80</td></tr> </tbody> </table> <p>The percentage ranges in the table below are based on the weight of the 4F Single Lumen (3.17g) and 6F Triple Lumen (7.26g) Pro-PICC® Valved Power Injectable PICCs.</p> <p style="text-align: center;">Pro-PICC® Power Injectable PICCs (Valved)</p> <table border="1"> <thead> <tr> <th>Material</th><th>% Weight (w/w)</th></tr> </thead> <tbody> <tr> <td>Polyurethane</td><td>86.58 - 91.51</td></tr> <tr> <td>Barium Sulfate</td><td>5.78 - 6.83</td></tr> <tr> <td>Silicone</td><td>2.64 - 6.83</td></tr> </tbody> </table> <p>Note: Accessories containing stainless steel may contain up to 0.4% weight of the CMR substance cobalt.</p>	Material	% Weight (w/w)	Polyurethane	58.88 - 64.09	Acetal Co-polymer	16.82 - 24.41	Acrylonitrile Butadiene Styrene	8.13 - 10.57	Barium Sulfate	2.82 - 11.80	Material	% Weight (w/w)	Polyurethane	86.58 - 91.51	Barium Sulfate	5.78 - 6.83	Silicone	2.64 - 6.83
Material	% Weight (w/w)																		
Polyurethane	58.88 - 64.09																		
Acetal Co-polymer	16.82 - 24.41																		
Acrylonitrile Butadiene Styrene	8.13 - 10.57																		
Barium Sulfate	2.82 - 11.80																		
Material	% Weight (w/w)																		
Polyurethane	86.58 - 91.51																		
Barium Sulfate	5.78 - 6.83																		
Silicone	2.64 - 6.83																		

	Note: The device should not be used if you are allergic to the above materials.	
Information on medicinal substances in the device	N/A	
How the device achieves its intended mode of action	<p>The subject devices utilize a Seldinger or Modified Seldinger technique to gain access. The main difference is one technique utilizes an Introduction Sheath and one does not. The Seldinger techniques for venous access are well-known surgical techniques used for inserting PICC devices. The instructions for use of each catheter are detailed in the IFUs. Catheters are to be inserted, manipulated and removed by a qualified, licensed physician or other qualified health care professional utilizing strict aseptic technique.</p> <p>Once in place, fluids are delivered or blood is withdrawn via the PICC catheter most commonly with a disposable tubing set or syringe. Catheter care includes use of a locking solution to maintain catheter patency. Catheter removal is normally done by gently pulling on the catheter, but removal may require that a surgical procedure be performed by a physician familiar with the appropriate techniques in some circumstances.</p>	
Sterilization Information	Contents sterile and non-pyrogenic in unopened, undamaged package. Sterilized by Ethylene Oxide.	
Description of accessories	Name of Accessory	Description of Accessory
	Guidewire	Acts as a path for other components.
	Introducer Needle	Placed into the target vein to gain access.
	Stylet	Assists in catheter placement.
	Peelable Introducer	Used to get central venous access.
	Scalpel	A cutting device.
	Syringe	Helps get blood return once the needle punctures the vein.

4. Risks and warnings

Contact your healthcare professional if you believe you are experiencing side-effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your healthcare professional if needed.

How potential risks have been controlled or managed	<p>There have been 670,138 devices sold since January 2019. There are side effects and risks associated with the device. These include:</p> <ul style="list-style-type: none"> • Infection • Bleeding • Device Removal • Device Replacement <p>These risks are reduced to an acceptable level. The labeling describes the risks. The benefit of the device is central venous</p>
---	--

	access when alternatives are not suitable. These benefits outweigh the risks.																																				
Remaining risks and undesirable effects	<p>The Pro-PICC® Power Injectable Peripherally Inserted Central Catheter is associated with risks. These include:</p> <ul style="list-style-type: none">• Procedural Delays• Thrombosis• Infections• Perforations• Embolism• Cardiac Event• Dissatisfaction <p>These risks are consistent with risks of other peripherally inserted central catheters. They are not unique to the Medcomp product. Some of the most common reactions include infection. Infection may be associated with general surgical procedure and hospitalization. Infection may not always be device-related.</p> <table><tr><th rowspan="4">Patient Residual Harm Category</th><th colspan="2">Quantification of Residual Risks</th></tr><tr><th>Complaints (01 January 2019 – 30 September 2024)</th><th>Post Market Clinical Follow-Up Activity Events</th></tr><tr><th>Units Sold: 670,138</th><th>Units Studied: 2,030</th></tr><tr><th># of Cases Per Event</th><th># of Cases Per Event</th></tr><tr><td>Allergic Reaction</td><td>Not Reported.</td><td>Not Reported.</td></tr><tr><td>Bleeding</td><td>1 Event in 650,000 Cases.</td><td>Not Reported.</td></tr><tr><td>Cardiac Event</td><td>Not Reported.</td><td>Not Reported.</td></tr><tr><td>Embolism</td><td>Not Reported.</td><td>Not Reported.</td></tr><tr><td>Infection</td><td>Not Reported.</td><td>1 Event in 50 Cases.</td></tr><tr><td>Perforation</td><td>Not Reported.</td><td>Not Reported.</td></tr><tr><td>Stenosis</td><td>Not Reported.</td><td>Not Reported.</td></tr><tr><td>Tissue Injury</td><td>1 Event in 650,000 Cases.</td><td>Not Reported.</td></tr><tr><td>Thrombosis</td><td>Not Reported.</td><td>Not Reported.</td></tr></table>	Patient Residual Harm Category	Quantification of Residual Risks		Complaints (01 January 2019 – 30 September 2024)	Post Market Clinical Follow-Up Activity Events	Units Sold: 670,138	Units Studied: 2,030	# of Cases Per Event	# of Cases Per Event	Allergic Reaction	Not Reported.	Not Reported.	Bleeding	1 Event in 650,000 Cases.	Not Reported.	Cardiac Event	Not Reported.	Not Reported.	Embolism	Not Reported.	Not Reported.	Infection	Not Reported.	1 Event in 50 Cases.	Perforation	Not Reported.	Not Reported.	Stenosis	Not Reported.	Not Reported.	Tissue Injury	1 Event in 650,000 Cases.	Not Reported.	Thrombosis	Not Reported.	Not Reported.
	Patient Residual Harm Category		Quantification of Residual Risks																																		
			Complaints (01 January 2019 – 30 September 2024)	Post Market Clinical Follow-Up Activity Events																																	
			Units Sold: 670,138	Units Studied: 2,030																																	
		# of Cases Per Event	# of Cases Per Event																																		
	Allergic Reaction	Not Reported.	Not Reported.																																		
	Bleeding	1 Event in 650,000 Cases.	Not Reported.																																		
	Cardiac Event	Not Reported.	Not Reported.																																		
	Embolism	Not Reported.	Not Reported.																																		
	Infection	Not Reported.	1 Event in 50 Cases.																																		
Perforation	Not Reported.	Not Reported.																																			
Stenosis	Not Reported.	Not Reported.																																			
Tissue Injury	1 Event in 650,000 Cases.	Not Reported.																																			
Thrombosis	Not Reported.	Not Reported.																																			
Warnings and precautions	<p>The below are warnings, precautions, or measures to be taken by patient:</p> <ul style="list-style-type: none">• Keep catheter dressing clean and dry. Ask your doctor for specific instructions on how to care for your catheter.• Avoid letting the catheter or catheter site go under water. Moisture near the catheter site can potentially lead to an infection. Patients must not swim, shower or soak dressing while bathing.• Contact your doctor right away if you notice any signs or symptoms of your catheter complications, such as																																				

	<ul style="list-style-type: none"> o The area around your line is getting red, swollen, bruised or warm to the touch. o Drainage from your catheter site. o The length of the catheter that sticks out of your insertion site gets longer. o Difficulty flushing your line because it seems to be blocked. • Avoid lifting heavy objects. • Do not have blood pressure readings taken on arm where the catheter is placed
Summary of any field safety correction action (FSCA)	There were no recalls for the device between 01 December 2023 to 30 September 2024.

5. Summary of clinical evaluation and post-market clinical follow-up

Clinical background of the device
<p>The subject devices have been available since 2007. The CE Mark was received in October 2007. US Food and Drug Administration (FDA) clearance was in September 2009. All models included are planned for distribution in the European Union.</p>
Clinical evidence for CE-marking
<p>The clinical literature review identified 16 articles relating to the safety and/or performance of the subject device when used as intended. These articles included approximately 4,047 cases. Two patient level data activities received information on 2,030 catheters. 37 user surveys have been received relating to this device.</p> <p>Findings from the clinical literature and data activities support the performance of the subject device. All data on the Pro-PICC® catheter has been evaluated. The benefits of the subject device outweigh the risks when the device is used as intended. The benefit of the device is facilitating blood withdraw for laboratory testing, delivery of fluids and medications for treatments including chemotherapy, and power injection of contrast media for CT examinations without requiring frequent needlesticks. These benefits are for patients in whom short-term or long-term peripheral access to the central venous system is deemed necessary based on the direction of a qualified, licensed physician.</p>
Safety
<p>There is sufficient data to prove conformity to the applicable requirements. The device is safe and performs as intended. The device is state of the art.</p> <p>Medcomp has reviewed:</p> <ul style="list-style-type: none"> • Post-Market Data • Medcomp Information Materials • Risk Management Documentation

The risks are appropriately displayed and consistent with the state of the art. The risks associated with the device are acceptable when weighed against the benefits.

There were 670,138 devices sold from January 1st 2019, to September 30th, 2024. Also, during this period there were 186 complaints received resulting in a 0.028% complaint frequency for the product family.

6. Possible therapeutic alternatives

When considering alternative treatments, it is recommended to contact your healthcare professional who can consider your individual situation. The Infusion Nurses Society (INS) Standards 2021 clinical practice guidelines have been used to support the below recommendations for treatments.

Therapy	Benefits	Disadvantages	Key Risks
Central Venous Catheters (CVCs)	<ul style="list-style-type: none"> • Easy access. • Minimizes repeat puncture. • Increased patient mobility. • Easier for outpatients. 	<ul style="list-style-type: none"> • Requires surgery. • Surgery risks. • Requires maintenance. • High risk of infection or thrombosis. 	<ul style="list-style-type: none"> • Infection • Occlusion • Malfunction • Thrombosis
Implantable Ports	<ul style="list-style-type: none"> • Less Vein Damage. • Easier to See and Access. • Reduces chance for corrosive medications to make skin contact • One puncture location. • Longer Dwell Time. • Can be permanent. 	<ul style="list-style-type: none"> • Requires surgery. • Surgery Risks. • Requires maintenance. 	<ul style="list-style-type: none"> • Infection • Embolism • Necrosis
Midline Catheters	<ul style="list-style-type: none"> • Patient comfort. • Longer dwell time than PIVs. • Lower risk of infection compared to IVs • No X-ray required. • Decreased chance of extravasation. 	<ul style="list-style-type: none"> • Not suitable for continuous injections of most vesicants or irritants 	<ul style="list-style-type: none"> • Phlebitis

Therapy	Benefits	Disadvantages	Key Risks
Peripherally Inserted Central Catheters (PICCs)	<ul style="list-style-type: none"> Decreased risk of catheter occlusion compared to CVC Fewer punctures compared to PIV 	<ul style="list-style-type: none"> Increased risk of deep vein thrombosis compared to CVC Pain/Discomfort over time Daily Life Adaption 	<ul style="list-style-type: none"> Deep vein thrombosis (DVT) Pulmonary embolism <ul style="list-style-type: none"> Venous thromboembolism (VTE) Post thrombotic syndrome
Peripheral Intravenous Catheters (PIVs)	<ul style="list-style-type: none"> No Surgery. 	<ul style="list-style-type: none"> Infection Bleeding Thrombosis Cannot be used for therapies with blistering agents Four days maximum use. 	<ul style="list-style-type: none"> Infection Phlebitis

7. Suggested training for users

The catheter should be inserted, manipulated, and removed by a qualified, licensed physician or other qualified health care professional under the direction of a physician.

Abbreviation	Definition
CE	Conformité Européenne (European Conformity)
cm	centimeter
CMR	Carcinogenic, mutagenic, reprotoxic
CT	Computerized Tomography (CAT Scan)
CVC	Central Venous Catheter
dba	Doing Business As
F	French (thickness of catheter)
FDA	Food and Drug Administration
FSCA	Field Safety Corrective Action
INS	Infusion Nurses Society
IV	Intravenous
N/A	Not Applicable
PA	Pennsylvania
PICC	Peripherally Inserted Central Catheter
PIV	Peripheral Intravenous Catheters
SSCP	Summary of Safety and Clinical Performance
USA	United States of America
w/w	Weight over Weight

Add copy to 'MDR Documentation' (Initial & Date): amm 24OCT2024