

Cath Only		
MR17035100	5F X 60CM SINGLE CT RATED PICC	5/BOX
MR17035200	5F X 55CM DUAL CT RATED PICC	5/BOX

Basic Set		
MR17033101	3F X 55CM SINGLE CT RATED PICC	5/BOX
MR17034101	4F X 55CM SINGLE CT RATED PICC	5/BOX
MR17035101	5F X 60CM SINGLE CT RATED PICC	5/BOX
MR17035201	5F X 55CM DUAL CT RATED PICC	5/BOX
MR17036201	6F X 60CM DUAL CT RATED PICC	5/BOX
MR17036301	6F X 60CM TRIPLE CT RATED PICC	5/BOX

Long Wire Set		
MR17033102	3F X 55CM SINGLE CT RATED PICC	5/BOX
MR17034102	4F X 55CM SINGLE CT RATED PICC	5/BOX
MR17035102	5F X 60CM SINGLE CT RATED PICC	5/BOX
MR17035202	5F X 55CM DUAL CT RATED PICC	5/BOX
MR17036202	6F X 60CM DUAL CT RATED PICC	5/BOX
MR17036302	6F X 60CM TRIPLE CT RATED PICC	5/BOX

Nursing Kit		
MR17033105	3F X 55CM SINGLE CT RATED PICC	5/BOX
MR17034105	4F X 55CM SINGLE CT RATED PICC	5/BOX
MR17035105	5F X 60CM SINGLE CT RATED PICC	5/BOX
MR17035205	5F X 55CM DUAL CT RATED PICC	5/BOX
MR17036205	6F X 60CM DUAL CT RATED PICC	5/BOX
MR17036305	6F X 60CM TRIPLE CT RATED PICC	5/BOX

Cut-Down Tray		
MRCTP41024	4F X 55CM SINGLE CT RATED PICC	5/BOX
MRCTP51024	5F X 60CM SINGLE CT RATED PICC	5/BOX
MRCTP52024	5F X 55CM DUAL CT RATED PICC	5/BOX
MRCTP62024	6F X 60CM DUAL CT RATED PICC	5/BOX

CATH ONLY INCLUDES

- (1) PICC w/ Stylet
- (1|2) Needle Free Valve(s)
- (1) Patient Information Pack
- (1) Patient Chart Sticker
- (1) Catheter Securement Device

BASIC SET INCLUDES

- (1) PICC w/ Stylet
- (1) Peelable Sheath Introducer
- (1) Mini Scalpel
- (1) Guidewire
- (1) 10cc Syringe
- (1) Tape Measure
- (1) Introducer Needle
- (1|2|3) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Patient Information Pack
- (1) Patient Chart Sticker

LONG WIRE SET INCLUDES

- (1) PICC
- (1) Peelable Sheath Introducer
- (1) Mini Scalpel
- (1) Guidewire
- (1) 10cc Syringe
- (1) Tape Measure
- (1) Introducer Needle
- (1|2|3) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Patient Information Pack
- (1) Patient Chart Sticker

NURSING KIT INCLUDES

- (1) PICC w/ Stylet
- (1) Peelable Sheath Introducer
- (1) Safety Scalpel
- (1) Guidewire
- (1) 10cc Syringe
- (1) Tape Measure
- (1) Safety Introducer Needle
- (1|2|3) Needle Free Valve(s)
- (1) Catheter Securement Device
- (1) Latex-Free Tourniquet
- (1) Patient Information Pack
- (1) Patient Chart Sticker

CUT-DOWN TRAY INCLUDES

- (1) PICC w/ Stylet
- (1) Adaptor w/ Sideport
- (1) Introducer Needle
- (1) Peelable Sheath Introducer
- (1) Safety Scalpel
- (1) Guidewire
- (1) Tape Measure
- (1) Catheter Securement Device

PRO-PICC^{CT}
CT RATED PICC CATHETER



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438 U.S.A.
Tel: 215-256-4201 | Fax: 215-256-1787
www.medcomp.net

CE2797

EC REP

ISO13485

Certified Quality System

MPS Medical Product Service GmbH
Borngasse 20
35619 Braunfels
Germany

PN2551 Rev. E 8/23

www.medcomp.net

PRO-PICC^{CT}

CT RATED PICC CATHETER



21 Gauge Introducer Needle with Echogenic Tip

- Super Sharp Needle




Guidewires

- Nitinol Material



Tearaway Sheath



TAPER

WHITE DEPTH MARKINGS

IMPORTANT RISK INFORMATION

Indications- for Use: The Pro-PICC[®] CT catheter is indicated for short or long term peripheral access to the central venous system for intravenous therapy and power injection of contrast media, and allows for central venous pressure monitoring. For blood sampling, infusion, or therapies, use a 4F or larger catheter. The maximum recommended infusion rate varies by catheter French size and is printed on the catheter.

Contraindications: The presence of device related infection, bacteremia, or septicemia is known or suspected. The patient's body size is insufficient to accommodate the size of the implanted device. The patient is known or is suspected to be allergic to materials contained in the device. There has been past irradiation of prospective insertion site. There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site. There are local tissue factors that may prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation