

Arch-flo™

medCOMP®

ORDERING INFORMATION

BASIC SET	CATALOG #	QTY/BOX
4F X 10CM ARCH-FLO™ CT MIDLINE, SINGLE LUMEN	MLCTAFD4S10	5

SETS CONTAIN:

- | | |
|------------------------------|--------------------------------|
| (1) Catheter | (1) Safety Scalpel |
| (1) Guidewire | (1) 10cc Syringe |
| (1) Needle Free Valve | (1) Tape Measure |
| (1) Peelable Introducer | (1) Catheter Securement Device |
| (1) Safety Introducer Needle | |



Arch-flo™


CT MIDLINE

Arch-flo™

Bridging the Gap

The Medcomp CT Midline Arch-Flo™ with SMARTip™ Technology provides short-term (1-4 weeks) peripheral access to the venous system for selected intravenous therapies and blood sampling. The use of Arch-Flo™ Bridges the Gap between PIVs and PICCs.

Our ENSUREID™ Catheter Identification System offers three points of reference to help clinicians identify the catheter as a midline. The unique combination of a yellow hub and luer system accented with bold “MIDLINE” font differentiates the device. Arch-Flo™ is placed using a familiar MST approach. Our tip design is what separates us from the standard tips.




- **Radius Tip** for a smooth surface.
- **Tip opening increased 150%** vs. standard tip design.

ENSUREID™



Optimized Taper



CT Rated
for 5cc/sec

MST
Placement