


DUO-FLOW®

SIDE X SIDE



FEATURES & BENEFITS

Double “D” internal lumen design

Tapered tip with large sideholes

Straight or Curved extensions

Available in 9F, 11F, 12F, & 14F to accommodate multiple patient populations

DUO-FLOW® SIDE X SIDE

Ordering Information

9F & 11F BASIC SET		
AS099101	9F X 10CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS099121	9F X 12CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS099151	9F X 15CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS099201	9F X 20CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AC099751	9F X 7.5CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC099121	9F X 12CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC099151	9F X 15CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AS0911121	11F X 12CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0911151	11F X 15CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0911201	11F X 20CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0911241	11F X 24CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AC0911101	11F X 10CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0911121	11F X 12CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0911151	11F X 15CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0911201	11F X 20CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX

9F/11F BASIC SET CONTAINS: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Dilator (1) Scalpel
(1) Adhesive Wound Dressing (2) End Caps

12F & 14F BASIC SET		
AS0912151	12F X 15CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0912201	12F X 20CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0912241	12F X 24CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AC0912131	12F X 13CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0912151	12F X 15CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0912201	12F X 20CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AS0914151	14F X 15CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0914201	14F X 20CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AS0914241	14F X 24CM DUO-FLOW® SIDE X SIDE, STRAIGHT	10/BOX
AC0914151	14F X 15CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0914201	14F X 20CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX
AC0914241	14F X 24CM DUO-FLOW® SIDE X SIDE, CURVED EXTENSION	10/BOX

12F/14F BASIC SET CONTAINS: (1) Catheter (1) Introducer Needle (1) Guidewire (2) Dilators (1) Scalpel
(1) Adhesive Wound Dressing (2) End Caps

IMPORTANT RISK INFORMATION

Indications for Use: The Duo-Flow® Side X Side double lumen catheter is intended for short-term central venous access for hemodialysis, apheresis and infusion.

Contraindications: The catheter is not intended for any purpose other than indicated in these instructions. The Duo-Flow® Side X Side double lumen catheter is intended for short-term (less than 30 days) use. Do not use this catheter in thrombosed vessels or for subclavian puncture when a ventilator is in use. Do not use this catheter when: The patient's body size is insufficient to accommodate the implanted device. The superficial or deep tissue will not permit adequate device stabilization and/or access. There are known physiological limitations that will not allow placement of the device. The patient has known or suspected allergies to any of the materials in the device. The patient has received significant radiation at the intended exit site or tunnel. The patient has severe chronic obstructive lung disease.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
215.256.4201
www.medcomp.net

PN2654 Rev. B 8/25