



DIGNITY® CT/ DIGNITY® MINI/ DIGNITY® DMP

Indications for Use: The CT Power Injectable Implantable Infusion Ports are indicated for patient therapies requiring repeated access to the vascular system. The port system can be used for infusion of medications, I.V. fluids, parenteral nutrition solutions, blood products, and for the withdrawal of blood samples. When used with a power injectable needle, the Power Injectable Implantable Infusion Port device is indicated for power injection of contrast media. For power injection of contrast media, the maximum recommended infusion rate is 5ml/s with a 19 or 20 gauge non-coring power injectable needle. The maximum recommended infusion rate is 2ml/s with a 22 gauge non-coring power injectable needle.

Contraindications: This device is contraindicated for catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates of pinch-off. The device is also contraindicated: When the presence of device related infection, bacteremia, or septicemia is known or suspected. When the patient's body size is insufficient for the size of the implanted device. When the patient is known or is suspected to be allergic to materials contained in the device. If severe chronic obstructive lung disease exists. If the prospective insertion site has been previously irradiated. If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures. If local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

Plastic
Port Material

**Pro-Fuse[®],
Pro-Fuse[®] LP,
Dignity[®] Mid
Sized,
Dignity[®] Mini,
Dignity[®] LP**
Configurations

Polyurethane
Catheter Material

**5cc/sec
@ 325psi**
Injection Rating

5F, 6.6F, 8F
French Sizes



21GA NEEDLE



PTFE SUPER SHEATH DILATOR



NITINOL MANDREL WIRE

PRO-FUSE [®] PORTS		
MRCTT80001DMP*	8F PRO-FUSE [®] DMP PORT KIT	5/BOX
MRCTT80041DMP	8F PRO-FUSE [®] LOW PROFILE DMP PORT KIT	5/BOX

DIGNITY [®] PORTS		
MRCTI50041MDMP	5F DIGNITY [®] MINI DMP PORT KIT	5/BOX
MRCTI50841MDMP	5F DIGNITY [®] MINI PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI50041DMP	5F DIGNITY [®] LOW PROFILE DMP PORT KIT	5/BOX
MRCTI50841DMP	5F DIGNITY [®] LOW PROFILE PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI50001DMP	5F DIGNITY [®] MID SIZED DMP PORT KIT	5/BOX
MRCTI50801DMP	5F DIGNITY [®] MID SIZED PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI66001DMP	6.6F DIGNITY [®] MID SIZED DMP PORT KIT	5/BOX
MRCTI66041MDMP	6.6F DIGNITY [®] MINI DMP PORT KIT	5/BOX
MRCTI66841MDMP	6.6F DIGNITY [®] MINI PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI66041DMP	6.6F DIGNITY [®] LOW PROFILE DMP PORT KIT	5/BOX
MRCTI66801DMP	6.6F DIGNITY [®] MID SIZED PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI66841DMP	6.6F DIGNITY [®] LOW PROFILE PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI80001DMP	8F DIGNITY [®] MID SIZED DMP PORT KIT	5/BOX
MRCTI80041MDMP	8F DIGNITY [®] MINI DMP PORT KIT	5/BOX
MRCTI80841MDMP	8F DIGNITY [®] MINI PRE ATTACHED DMP PORT KIT	5/BOX
MRCTI80041DMP	8F DIGNITY [®] LOW PROFILE DMP PORT KIT	5/BOX
MRCTI80841DMP	8F DIGNITY [®] LOW PROFILE PRE-ATTACHED DMP PORT KIT	5/BOX
MRCTI80801DMP	8F DIGNITY [®] MID SIZED PRE-ATTACHED DMP PORT KIT	5/BOX

- KITS CONTAIN:
- (1) Catheter
 - (1) Peelable Introducer
 - (1) Scalpel
 - (1) 21ga Needle
 - (1) .018" Nitinol Guidewire
 - (1) Blunt Tip Needle
 - (1) Tunneler
 - (1) CT Implantable Port
 - (1) 22ga Huber Needle
 - (1) 10cc Syringe
 - (2) Catheter Locks
 - (1) Patient Chart Sticker

*MRCTT80001DMP has Silicone Filled Suture Holes

Refer to the Table of Contents for Important Risk Information regarding this device.