



Power Injectable Peripherally Inserted Central Catheter

Basic Set		
MR17033101	3F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17034101	4F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035101	5F X 60CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035201	5F X 55CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036201	6F X 60CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036301	6F X 60CM TRIPLE LUMEN PRO-PICC®	5/BOX

Longwire Set		
MR17033102	3F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17034102	4F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035102	5F X 60CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035202	5F X 55CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036202	6F X 60CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036302	6F X 60CM TRIPLE LUMEN PRO-PICC®	5/BOX

Nursing Set		
MR17033105	3F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17034105	4F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035105	5F X 60CM SINGLE LUMEN PRO-PICC®	5/BOX
MR17035205	5F X 55CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036205	6 X 60CM DUAL LUMEN PRO-PICC®	5/BOX
MR17036305	6 X 60CM TRIPLE LUMEN PRO-PICC®	5/BOX

Full Cut Down Set		
MRCTP31024	3F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MRCTP41024	4F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MRCTP51024	5F X 60CM SINGLE LUMEN PRO-PICC®	5/BOX
MRCTP52024	5F X 55CM DUAL LUMEN PRO-PICC®	5/BOX
MRCTP62024	6F X 60CM DUAL LUMEN PRO-PICC®	5/BOX
MRCTP63024	6F X 60CM TRIPLE LUMEN PRO-PICC®	5/BOX

Full Cut Down Set with 70cm Guidewire		
MRCTP41028	4F X 55CM SINGLE LUMEN PRO-PICC®	5/BOX
MRCTP51028	5F X 60CM SINGLE LUMEN PRO-PICC®	5/BOX
MRCTP52028	5F X 55CM DUAL LUMEN PRO-PICC®	5/BOX
MRCTP62028	6F X 60CM DUAL LUMEN PRO-PICC®	5/BOX

*3F Kit has a 3.5F Peelable Introducer
4F Kit has a 4.5F Peelable Introducer
5F Kit has a 5.5F Peelable Introducer
6F Kit has a 6.5F Peelable Introducer

Basic Set Contains:

- (1) Catheter w/ Stylet
- (1) Adaptor w/ Sideport
- (1) Peelable Introducer*
- (1) Introducer Needle
- (1) Syringe
- (1) Guidewire
- (1 | 2 | 3) Needleless Connector(s)
- (1) Securement Device
- (1) Scalpel
- (1) Tape Measure
- (1) Patient ID Card
- (1) Patient Information Packet

Longwire Set Contains:

- (1) Catheter
- (1) Peelable Introducer*
- (1) Introducer Needle
- (1) Syringe
- (1) Guidewire
- (1 | 2 | 3) Needleless Connector(s)
- (1) Securement Device
- (1) Scalpel
- (1) Tape Measure
- (1) Patient ID Card
- (1) Patient Information Packet

Nursing Set Contains:

- (1) Catheter w/ Stylet
- (1) Adaptor w/ Sideport
- (1) Tourniquet
- (1) Peelable Introducer*
- (1) Introducer Needle
- (1) Syringe
- (1) Guidewire
- (1 | 2 | 3) Needleless Connector(s)
- (1) Securement Device
- (1) Scalpel
- (1) Tape Measure
- (1) Patient ID Card
- (1) Patient Information Packet

Full Cut Down Set Contains:

- (1) Catheter w/ Stylet
- (1) Adaptor w/ Sideport
- (1) Peelable Introducer*
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Securement Device
- (1) Tape Measure
- (1) Patient ID Card
- (1) Patient Information Packet

Full Cut Down Set with 70cm Guidewire Contains:

- (1) Catheter w/ Stylet
- (1) Adaptor w/ Sideport
- (1) Peelable Introducer*
- (1) Introducer Needle
- (1) Guidewire
- (1) Scalpel
- (1) Securement Device
- (1) Tape Measure
- (1) Patient ID Card
- (1) Patient Information Packet



Power Injectable Peripherally Inserted Central Catheter



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201
www.medcomp.net

ISO13485
Certified Quality System

PN2690 Rev. A 9/26



Power Injectable Peripherally Inserted Central Catheter

 **medCOMP**
www.medcomp.net

Single, Dual, and Triple Lumen Options




21 Gauge Introducer Needle with Echogenic Tip
• Super Sharp Needle



Guidewires
• Nitinol Material



Tearaway Sheath



IMPORTANT RISK INFORMATION

Indications For Use: The Pro-PICC® Power Injectable Peripherally Inserted Central Catheter is indicated for short-term or long-term peripheral access to the central venous system for blood sampling, intravenous administration of fluids or medications, central venous pressure monitoring and power injection of contrast media.

Contraindications:

- The presence of device related local infection, bacteremia or septicemia is known or suspected.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- There are local tissue factors that may prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation