

**Pro-PICC® Valved Power Injectable
Peripherally Inserted Central Catheter**

CODE	DESCRIPTION	UOM
MR82034101	4F X 55CM SINGLE LUMEN VALVED CT PICC – BASIC SET	5/BOX
MR82035101	5F X 60CM SINGLE LUMEN VALVED CT PICC – BASIC SET	5/BOX
MR82035201	5F X 55CM DUAL LUMEN VALVED CT PICC – BASIC SET	5/BOX
MR82036301	6F X 60CM TRIPLE LUMEN VALVED CT PICC – BASIC SET	5/BOX

Basic Set Contains: (1) Catheter (1) Peelable Introducer* (1) Introducer Needle (1) Syringe (1) Guidewire (1) Needleless Connector(s) (1) Securement Device (1) Scalpel (1) Tape Measure (1) Patient Information Packet (1) Patient ID Card

*4F Kit has a 4.5F Peelable Introducer, 5F Kit has a 5.5F Peelable Introducer, 6F Kit has a 6.5F Peelable Introducer.

Important Risk Information

Indications For Use: The Pro-PICC® Valved Power Injectable Peripherally Inserted Central Catheter is indicated for short-term or long-term peripheral access to the central venous system for blood sampling, intravenous administration of fluids or medications, central venous pressure monitoring and power injection of contrast media.


Contraindications:

- The presence of device related infection, bacteremia, or septicemia is known or suspected.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- There has been past irradiation of prospective insertion site.
- There have been previous episodes of venous thrombosis or vascular surgical procedures at the prospective placement site.
- There are local tissue factors that may prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

PRO-PICC®
VALVED POWER INJECTABLE
PERIPHERALLY INSERTED
CENTRAL CATHETERS





Features & Benefits

- **Valve Identification**
All valved PICCs are labeled on the extensions.
- **Power Injection Rates**
5cc/sec at 325psi on all Pro-PICC® valved catheters.
- **Saline Only Flush**
- **Clamp Free Extensions**

