

CT MIDLINE

Power Injectable Midline

FEATURES & BENEFITS

TIP DESIGN

Trimable

AVAILABLE IN

3F, 4F & 5F

POWER INJECTION RATING

3F Single- 3cc/sec

4F & 5F Single- 5cc/sec

5F Dual- 7cc/sec

LENGTH

15cm, 20cm

CONFIGURATION

Single & Dual




Markings for catheter trimming



Printed identification denotes maximum flow rates for contrast injection.

With MST Placement



CT MIDLINE


Power Injectable Midline

BASIC SET

MLCT3S15	3F X 15CM SINGLE CT MIDLINE	5/BOX
MLCT4S20	4F X 20CM SINGLE CT MIDLINE	5/BOX
MLCT5S20	5F X 20CM SINGLE CT MIDLINE	5/BOX
MLCT5D20	5F X 20CM DUAL CT MIDLINE	5/BOX

Basic Set Contains: (1) Catheter With Stylet (1) Adaptor With Sideport (1) Peelable Introducer* (1) Introducer Needle (1) Syringe (1) Guidewire (1) 2 Needleless Connector(s) (1) Securement Device (1) Scalpel (1) Tape Measure (1) Patient Information Packet (1) Patient ID Card

*3F Kit has a 3.5F Peelable Introducer
4F Kit has a 4.5F Peelable Introducer
5F Kit has a 5.5F Peelable Introducer



IMPORTANT RISK INFORMATION

Indications for Use: The CT Midline Power Injectable Midline is indicated for short-term peripheral access to the peripheral venous system for intravenous administration of fluids or medications and power injection of contrast media.

Contraindications:

- This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels.
- The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.)
- The presence of device related bacteremia or septicemia.
- History of mastectomy on insertion side.
- Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site.
- Fever of unknown origin.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Past irradiation of prospective insertion site.
- Local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.

1499 Delp Drive, Harleysville, PA 19438, U.S.A.

P: 215.256.4201

www.medcomp.net

ISO13485

Certified Quality System

PN2700 Rev. A 9/26