




Important Risk Information

Indications For Use: The 1.9F and 2.6F VASCU-PICC® Peripherally Inserted Central Catheter is indicated for short-term or long-term peripheral access to the central venous system for the intravenous administration of fluids or medications.

Contraindications:

- This catheter is not intended for any use other than that, which is indicated. Do not implant catheter in thrombosed vessels.
- The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.)
- The presence of device related bacteremia or septicemia.
- Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site.
- Fever of unknown origin.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Past irradiation of prospective insertion site.
- Local tissue factors will prevent proper device stabilization and/or access.
- Known tape or Zinc Oxide adhesive allergies.
- This catheter is not suitable for insertion through non-superficial veins.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation


1.9F & 2.6F VASCU-PICC®

Peripherally Inserted Central Catheter



1.9F & 2.6F VASCU-PICC®

Peripherally Inserted Central Catheter



Basic Set With Cuff		
MR17012608	2.6F X 50CM DUAL LUMEN	5/BOX

Basic Set With Cuff Contains:
(1) Catheter With Cuff and Stylet (1) Adaptor With Sideport (1) OTN Tearaway Introducer* (10) Gauze (1) Syringe (2) Needleless Connectors (1) Tourniquet (1) Tape Measure (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

Catheter Set		
MR17011100	1.9F X 50CM SINGLE LUMEN	10/BOX
MR17012600	2.6F X 50CM DUAL LUMEN	10/BOX

Catheter Set Contains:
(1) Catheter With Stylet (1) Adaptor With Sideport (1 | 2) Needleless Connector(s) (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

Basic Set		
MR17011101	1.9F X 50CM SINGLE LUMEN	5/BOX
MR17011121	1.9F X 20CM SINGLE LUMEN	5/BOX
MR17012601	2.6F X 50CM DUAL LUMEN	5/BOX
MR17012621	2.6F X 20CM DUAL LUMEN	5/BOX

Basic Set Contains:
(1) Catheter With Stylet (1) Adaptor With Sideport (1) OTN Tearaway Introducer* (1 | 2) Needleless Connector(s) (10) Gauze (1) Syringe (1) Tourniquet (1) Tape Measure (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

Basic MST Set		
MR170111024S	1.9F X 50CM SINGLE LUMEN	5/BOX
MR170111224S	1.9F X 20CM SINGLE LUMEN	5/BOX
MR170126024S	2.6F X 50CM DUAL LUMEN	5/BOX
MR170126224S	2.6F X 20CM DUAL LUMEN	5/BOX

Basic MST Set Contains:
(1) Catheter With Stylet (1) Guidewire (1) Adaptor With Sideport (1) Peelable Introducer** (1 | 2) Needleless Connector(s) (1) IV Catheter (10) Gauze (1) Syringe (1) Tourniquet (1) Tape Measure (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

Basic Set Without Stylet		
VP1.9S20-NS	1.9F X 20CM SINGLE LUMEN	5/BOX
VP1.9S50-NS	1.9F X 50CM SINGLE LUMEN	5/BOX

Basic Set Without Stylet Contains:
(1) Catheter (1) OTN Tearaway Introducer* (1) Needleless Connector (10) Gauze (1) Syringe (1) Tourniquet (1) Tape Measure (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

*1.9F Kit has a 2F OTN Tearaway Introducer
2.6F Kit has a 3F OTN Tearaway Introducer

**1.9F Kit has a 2F Peelable Introducer
2.6F Kit has a 3F Peelable Introducer