

IMPORTANT RISK INFORMATION

Indications For Use: The Vascu-PICC® Peripherally Inserted Central Catheter Instructions are designed for Short or Long Term central venous catheterization (intravenous administration of fluids, medications, and/or when nutritional therapy is prescribed).

Contraindications:

- This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels.
- The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.)
- The presence of device related bacteremia or septicemia.
- History of mastectomy on insertion side.
- Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site.
- Fever of unknown origin.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Past irradiation of prospective insertion site.
- Local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



VASCU-PICC®
PERIPHERALLY INSERTED CENTRAL CATHETERS



VASCU-PICC®

PERIPHERALLY INSERTED CENTRAL CATHETERS

Ranging in size from 3F to 6F, and available with tapered or non-tapered lumens, Medcomp's Vascu-PICC® line of catheters provides a full range of options designed to meet the unique access challenges of patients young and old.



ORDERING INFORMATION

VASCU-PICC® BASIC SET		
MR17013101	3F X 60CM SINGLE LUMEN	5/BOX
MR17014101	4F X 60CM SINGLE LUMEN	5/BOX
MR17015101	5F X 60CM SINGLE LUMEN	5/BOX
MR17014201	4F X 60CM DUAL LUMEN	5/BOX
MR17015201	5F X 60CM DUAL LUMEN	5/BOX
MR17016201	6F X 60CM DUAL LUMEN	5/BOX

Set Contains: (1) Catheter w/ Stylet (1) Adaptor w/ Sideport (1) Peelable Introducer* (1) Needleless Connector(s) (1) Introducer Needle (1) Syringe (1) Guidewire (1) Securement Device (1) Safety Scalpel (1) Tape Measure (1) Patient Information Packet (1) Patient ID Card

VASCU-PICC® NURSING SET		
MR17013105	3F X 60CM SINGLE LUMEN	5/BOX
MR17014105	4F X 60CM SINGLE LUMEN	5/BOX
MR17015105	5F X 60CM SINGLE LUMEN	5/BOX
MR17014205	4F X 60CM DUAL LUMEN	5/BOX
MR17015205	5F X 60CM DUAL LUMEN	5/BOX

Set Contains: (1) Catheter w/ Stylet (1) Adaptor w/ Sideport (1) Tourniquet (1) Peelable Introducer* (1) Needleless Connector(s) (1) Safety Introducer Needle (1) Guidewire (1) Syringe (1) Safety Scalpel (1) Securement Device (1) Tape Measure (1) Patient Information Packet (1) Patient ID Card

*3F Kit has a 3.5F Peelable Introducer, 4F Kit has a 4.5F Peelable Introducer, 5F Kit has a 5.5F Peelable Introducer, 6F Kit has a 6.5F Peelable Introducer

