




FEATURES & BENEFITS

LUMENS
Single & Dual

AVAILABLE IN
3F, 4F & 5F

Cuff Placement at: 2cm or 5cm.





5CM FROM HUB-BASIC SET

MR28013151	3F X 60CM SINGLE VASCU-LINE® WITH CUFF	5/BOX
MR28014151	4F X 60CM SINGLE VASCU-LINE® WITH CUFF	5/BOX
MR28014251	4F X 60CM DUAL VASCU-LINE® WITH CUFF	5/BOX
MR28015251	5F X 60CM DUAL VASCU-LINE® WITH CUFF	5/BOX

2CM FROM HUB-BASIC SET

MR28013121	3F X 60CM SINGLE VASCU-LINE® WITH CUFF	5/BOX
MR28014221	4F X 60CM DUAL VASCU-LINE® WITH CUFF	5/BOX

Basic Set Contains: (1) Catheter with Stylet (1) Peelable Introducer* (1) Introducer Needle (1) Syringe (1) Guidewire (1) Needleless Connector(s) (1) Securement Device (1) Scalpel (1) Tape Measure (1) Tunneler (1) Patient Information Packet (1) Patient ID Card

*3F kit has a 3.5F Peelable Introducer
 4F kit has a 4.5F Peelable Introducer
 5F kit has a 5.5F Peelable Introducer

IMPORTANT RISK INFORMATION

Indications for Use:

The Vascu-Line® Central Venous Catheter is indicated for short-term or long-term access to the central venous system for intravenous administration of fluids or medications.

Contraindications:

- The presence of device related infection, bacteremia, or septicemia is known or suspected.
- This catheter is intended for short or long-term vascular access and should not be used for any purpose other than indicated in these instructions.
- The patient is known or is suspected to be allergic to materials contained in the device.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.
 1499 Delp Drive, Harleysville, PA 19438, U.S.A.
 P: 215.256.4201
www.medcomp.net

ISO13485
 Certified Quality System

PN2705 Rev. A 9/26