

# HEMO-FLOW<sup>®</sup>

LONG-TERM HEMODIALYSIS CATHETER



## Features & Benefits

- › Double “D” Internal Lumen Design
  - › Step Tip Design
- › Straight or Pre-Curved Lumen Configurations
  - › Available in 14.5F

**14.5F HEMO-FLOW®**  
LONG-TERM HEMODIALYSIS CATHETER

CATHETER SET		
HFS24E.	14.5F X 24CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS28E.	14.5F X 28CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS32E.	14.5F X 32CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS36E.	14.5F X 36CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS40E.	14.5F X 40CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS 55	14.5F X 55CM HEMO-FLOW® CATHETER, STRAIGHT	5/BOX
HFS24PCE.	14.5F X 24CM HEMO-FLOW® CATHETER, PRE-CURVED	5/BOX
HFS28PCE.	14.5F X 28CM HEMO-FLOW® CATHETER, PRE-CURVED	5/BOX
HFS32PCE.	14.5F X 32CM HEMO-FLOW® CATHETER, PRE-CURVED	5/BOX
HFS36PCE.	14.5F X 36CM HEMO-FLOW® CATHETER, PRE-CURVED	5/BOX

**Set Contents:** (1) Catheter (1) Introducer Needle (1) Guidewire w/Advancer (1) Tunneler (1) Dilator (1) Valved Peelable Introducer (1) Scalpel (2) End Cap (1) Patient ID Card (1) Patient Information Packet

**Indications for Use:** The Hemo-Flow® Catheter are indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

**Contraindications:** Known or suspected allergies to any of the components of the catheter or kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit [www.medcomp.net/patientinformation](http://www.medcomp.net/patientinformation)

