



Power Injectable Central Venous Catheter

FEATURES & BENEFITS

LUMENS

Single, Dual & Triple

POWER INJECTION RATING

5cc/sec

AVAILABLE IN

5F, 6F & 7F



Cuff Placement at:
2cm or 5cm.



Extension legs are
marked Power Injectable.





Power Injectable Central Venous Catheter

BASIC SET (CUFF 5CM FROM HUB)		
MR28035101	5F X 60CM SINGLE PRO-LINE® CT	5/BOX
MR28036101	6F X 60CM SINGLE PRO-LINE® CT	5/BOX
MR28037101	7F X 60CM SINGLE PRO-LINE® CT	5/BOX
MR28035201	5F X 55CM DUAL PRO-LINE® CT	5/BOX
MR28036201	6F X 60CM DUAL PRO-LIN® CT	5/BOX
MR28037201	7F X 60CM DUAL PRO-LINE® CT	5/BOX
MR28036301	6F X 60CM TRIPLE PRO-LINE® CT	5/BOX
BASIC SET (CUFF 2CM FROM HUB)		
MR28035121	5F X 60CM SINGLE PRO-LINE® CT	5/BOX
MR28036121	6F X 60CM SINGLE PRO-LINE® CT	5/BOX
MR28035221	5F X 55CM DUAL PRO-LINE® CT	5/BOX
MR28036221	6F X 60CM DUAL PRO-LINE® CT	5/BOX


Basic Set Contains: (1) Catheter With Stylet (1) Introducer Needle (1) Peelable Introducer* (1) Syringe (1) Guidewire (1) Scalpel (1) Tape Measure (2) Tunnelers (1 I 2 I 3) Needleless Connector(s) (1) Securement Device (1) Patient Information Packet (1) Patient ID Card

*5F kit has a 5.5F Peelable Introducer

6F Single Lumen Kit has a 6F Peelable Introducer

6F Dual & 6F Triple Lumen Kit has a 6.5F Peelable Introducer

7F Kit has a 7F Peelable Introducer



IMPORTANT RISK INFORMATION

Indications for Use:

The Pro-Line® Power Injectable Central Venous Catheter is indicated for short-term or long-term access to the central venous system for intravenous administration of fluids or medications and power injection of contrast media.

Contraindications:

- The presence of device related infection, bacteremia, or septicemia is known or suspected.
- This catheter is intended for short or long-term vascular access and should not be used for any purpose other than indicated in these instructions.
- The patient is known or is suspected to be allergic to materials contained in the device.

Refer to Instructions for Use for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.

1499 Delp Drive, Harleysville, PA 19438, U.S.A.

P: 215.256.4201

www.medcomp.net

ISO13485
Certified Quality System

PN2715 Rev. A 9/26