



DUO-FLOW[®] 400XL


SHORT-TERM HEMODIALYSIS CATHETER




FEATURES & BENEFITS




INTERNAL LUMEN DESIGN
Coaxial



MATERIAL
Polyurethane

TIP DESIGN
Tapered


FRENCH SIZE
14F


CONFIGURATION

Straight 	Pre-Curved 	Curved Extensions 
---	--	--



DUO-FLOW[®] 400XL

SHORT-TERM HEMODIALYSIS CATHETER

BASIC SET		
DFXL144MT	14F X 12CM DUO-FLOW [®] 400XL, STRAIGHT	10/BOX
DFXL146MTE	14F X 15CM DUO-FLOW [®] 400XL, STRAIGHT	10/BOX
DFXL148MTE	14F X 20CM DUO-FLOW [®] 400XL, STRAIGHT	10/BOX
DFXL149MTE	14F X 24CM DUO-FLOW [®] 400XL, STRAIGHT	10/BOX
DFXL144IJSE	14F X 12CM DUO-FLOW [®] 400XL, PRE-CURVED	10/BOX
DFXL146IJSE	14F X 15CM DUO-FLOW [®] 400XL, PRE-CURVED	10/BOX
DFXL148IJSE	14F X 20CM DUO-FLOW [®] 400XL, PRE-CURVED	10/BOX
DFXL149IJSE	14F X 24CM DUO-FLOW [®] 400XL, PRE-CURVED	10/BOX
DFXL148CES	14F X 20CM DUO-FLOW [®] 400XL, CURVED EXTENSIONS	10/BOX

Basic Set Contents: (1) Catheter (1) Guidewire (1) Introducer Needle (1 | 2*) Dilator (2) End Cap (1) Adhesive Wound Dressing*

*DFXL144MT & DFXL148CES Only

IMPORTANT RISK INFORMATION

Indications for Use: The Medcomp[®] Duo-Flow and Duo-Flow 400XL Double Lumen Catheters are designed for acute hemodialysis and apheresis. They may be inserted percutaneously and are ideally placed in the internal jugular vein. Although these catheters may be inserted into the subclavian or femoral vein, the internal jugular is the preferred site.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



Medical Components, Inc.
 1499 Delp Drive, Harleysville, PA 19438, U.S.A.
 P: 215.256.4201 F: 215.256.1787
www.medcomp.net

ISO13485
 Certified Quality System

PN2725BSI Rev. A 5/24