

FREE FLOW ST

SHORT-TERM HEMODIALYSIS CATHETER



FEATURES & BENEFITS

INTERNAL LUMEN DESIGN

Coaxial



MATERIAL

Polyurethane

TIP DESIGN

Coaxial



FRENCH SIZE

12.5F, 15.5F




CONFIGURATIONS

Straight



IJ



FREE FLOW ST

SHORT-TERM HEMODIALYSIS CATHETER

12.5F BASIC SET		
MFFS1215S	12.5F X 15CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1220S	12.5F X 20CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1224S	12.5F X 24CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1212IJ	12.5F X 12CM FREE FLOW ST, IJ	5/BOX
MFFS1215IJ	12.5F X 15CM FREE FLOW ST, IJ	5/BOX
MFFS1220IJ	12.5F X 20CM FREE FLOW ST, IJ	5/BOX
11F BASIC SET		
MFFS1515S	15.5F X 15CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1520S	15.5F X 20CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1524S	15.5F X 24CM FREE FLOW ST, STRAIGHT	5/BOX
MFFS1512IJ	15.5F X 12CM FREE FLOW ST, IJ	5/BOX
MFFS1515IJ	15.5F X 15CM FREE FLOW ST, IJ	5/BOX
MFFS1520IJ	15.5F X 20CM FREE FLOW ST, IJ	5/BOX

Basic Sets Contents: (1) Catheter with Stylet (1) Guidewire (1) Guidewire Advancer
(1) Introducer Needle (3) Dilator (2) End Cap

IMPORTANT RISK INFORMATION

Indications for Use: The ST Catheter is indicated for use in attaining Short-Term vascular access for Hemodialysis, Hemofiltration, and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites include the subclavian and femoral veins as required. This catheter is indicated for a duration less than (30) days. For femoral placement, monitor catheter condition closely.

Contraindications: This catheter is intended for Short-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

