

# PERITONEAL DIALYSIS CATHETERS

## Accessories

MBCA




MTC



MCS



MTS



### PD ACCESSORIES

MBCA	BARBED LUER LOCK ADAPTOR	20/BOX
MTC	TITANIUM BARBED LUER LOCK ADAPTOR	5/BOX
MCS	PD CATHETER STYLET	5/BOX
MTS	PD CATHETER TUNNELER	10/BOX
2758BSI	PD STENCIL MAP PACKET	EACH

### IMPORTANT RISK INFORMATION

#### I-SERIES & V-SERIES

**Indications for Use:** The I-Series / V-Series Peritoneal Dialysis Catheter is indicated for short-term or long-term use where peritoneal access is required for the purpose of peritoneal dialysis for Acute Kidney Injury (AKI) or Chronic Kidney Disease (CKD).

**Contraindications:** The catheter is not intended for the Tenckhoff trocar method of insertion. The catheter is contraindicated for patients with: Infected anterior abdominal wall. Active peritonitis. Known or suspected allergies to any of the components of the catheter or the kit.

#### BARBED LUER LOCK ADAPTOR

**Indications for Use:** The Medcomp® Barbed Luer Lock Adaptor is indicated for long-term use as a connection point for peritoneal dialysis catheters and male luers.

#### TITANIUM LUER ADAPTOR

**Indications for Use:** The Medcomp® Titanium Barbed Luer Lock Adaptor is indicated for long-term use as a connection point for peritoneal dialysis catheters and male luers.

#### PD CATHETER STYLET

**Indications for Use:** The Medcomp® Peritoneal Catheter Stylet is indicated for transient use to keep a peritoneal dialysis catheter straight and maneuverable during insertion into the peritoneum.

#### PD CATHETER TUNNELER

**Indications for Use:** The Medcomp® Peritoneal Catheter Tunneler is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening for the placement of a Medcomp® peritoneal dialysis catheter.

#### PD ACCESSORY CONTRAINDICATIONS


Known or suspected allergies to the device.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit [www.medcomp.net/patientinformation](http://www.medcomp.net/patientinformation)



V-Series



I-Series



Medical Components, Inc.  
1499 Delp Drive, Harleysville, PA 19438, U.S.A.  
P: 215.256.4201 F: 215.256.1787  
[www.medcomp.net](http://www.medcomp.net)

ISO13485  
Certified Quality System

PN2744BSI REV. A 5/24

medCOMP

[www.medcomp.net](http://www.medcomp.net)

## I-Series

Material  
**Silicone**

**SIL**

Tip Design  
**Straight & Coiled**



French Size  
**15F**



Cuffs  
**Single & Dual**



### 15F CATHETER ONLY STRAIGHT

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)
MPP41S	15F X 41CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	15	-
MSC31S	15F X 31CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5/BOX	10.25	-
MSC37S	15F X 37CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5/BOX	16.25	-
MSC42S	15F X 42CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5/BOX	21	-
MSC46S	15F X 46CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5/BOX	25	-
MDC31S	15F X 31CM I-SERIES CATHETER, DUAL CUFF	5/BOX	6.5	3
MDC32S	15F X 32.25CM I-SERIES CATHETER, DUAL CUFF	5/BOX	7.25	3.5
MDC37S	15F X 37CM I-SERIES CATHETER, DUAL CUFF	5/BOX	10.5	5
MDC42S	15F X 42CM I-SERIES CATHETER, DUAL CUFF	5/BOX	15	5
MDC47S	15F X 47CM I-SERIES CATHETER, DUAL CUFF	5/BOX	20	5

### 15F CATHETER ONLY COILED

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)	Cuff to Coil (CM)
MPP40C	15F X 40.25CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	26	-	16.36
MPP57C	15F X 57.5CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	31	-	13
MPP60C	15F X 60CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	34	-	24.45
MDC57C	15F X 57.5CM I-SERIES CATHETER, DUAL CUFF	5/BOX	31	6	12
MDC63C	15F X 63CM I-SERIES CATHETER, DUAL CUFF	5/BOX	38.25	6	17.68

### 15F BASIC SET STRAIGHT

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)
MPD-118	15F X 18CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	8	-
MPD-141	15F X 41CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	15	-
MPD-146	15F X 46CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5/BOX	20	-
MPD-146S	15F X 46CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5/BOX	25	-
MPD-237	15F X 37CM I-SERIES CATHETER, DUAL CUFF	5/BOX	10.5	5
MPD-242	15F X 42CM I-SERIES CATHETER, DUAL CUFF	5/BOX	15	5

### 15F BASIC SET COILED

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)	Cuff to Coil (CM)
MPD-160	15F X 60CM I-SERIES, SINGLE PREPERITONEAL CUFF	5/BOX	34	-	24.45
MPD-257	15F X 57.5CM I-SERIES, DUAL CUFF	5/BOX	31	6	12
MPD-263	15F X 63CM I-SERIES, DUAL CUFF	5/BOX	38.25	6	17.68

**Catheter Only Set Contents:** (1) Catheter (1) Female Luer (1) Clamp (1) End Cap

**Basic Set Contents:** (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (3) Gauze (1) End Cap

## V-Series

Material  
**Silicone**

**SIL**

Tip Design  
**Straight & Coiled**



French Size  
**15F**



Cuffs  
**Dual**



### 15F CATHETER ONLY STRAIGHT

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)
MC200301	15F X 43CM V-SERIES CATHETER, LEFT	5/BOX	14.5	5
MC200302	15F X 43CM V-SERIES CATHETER, RIGHT	5/BOX	14.5	5

### 15F CATHETER ONLY COILED

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)	Cuff to Coil (CM)
MC200303	15F X 39CM V-SERIES CATHETER, LEFT	5/BOX	16.5	3	6.94
MC200401	15F X 63CM V-SERIES CATHETER, LEFT	5/BOX	34	5	17.82
MC200402	15F X 63CM V-SERIES CATHETER, RIGHT	5/BOX	34	5	17.82

### 15F BASIC SET STRAIGHT

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)
MC200101	15F X 43CM V-SERIES CATHETER, LEFT	5/BOX	14.5	5
MC200102	15F X 43CM V-SERIES CATHETER, RIGHT	5/BOX	14.5	5

### 15F BASIC SET COILED

Code	Description	Box Qty	Distal Tip to Cuff (CM)	Cuff to Cuff (CM)	Cuff to Coil (CM)
MC200201	15F X 63CM V-SERIES CATHETER, LEFT	5/BOX	34	5	17.82
MC200202	15F X 63CM V-SERIES CATHETER, RIGHT	5/BOX	34	5	17.82

**Catheter Only Set Contents:** (1) Catheter (1) Female Luer (1) Clamp (1) End Cap

**Basic Set Contents:** (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (3) Gauze (1) End Cap