


PERITONEAL DIALYSIS CATHETERS

Accessories



PD Accessories		
MBCAGL	BARBED LUER LOCK ADAPTOR	20/BOX
MTCGL	TITANIUM BARBED LUER LOCK ADAPTOR	5/BOX
TBCGL	TITANIUM BARBED CONNECTOR	50/BOX
MCSGL	PD CATHETER STYLET	5/BOX
MTSGL	PD CATHETER TUNNELER	10/BOX
FSXGL	X-SERIES CURVED TUNNELER	5/BOX
TRXGL	X-SERIES TUNNELING ROD	5/BOX

IMPORTANT RISK INFORMATION

I-SERIES & V-SERIES

Indications for Use: The Peritoneal Catheter is indicated for chronic peritoneal dialysis.

Contraindications: This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy.

X-SERIES

Indications for Use: The Peritoneal Catheter is indicated for chronic peritoneal dialysis. The extended length of the X-Series Catheter makes it especially applicable for peritoneal dialysis patients when it is necessary to locate the skin exit site remote from the usual lower abdominal region. The catheter may be particularly indicated in patients with obesity, floppy abdominal skin folds, urinary or fecal incontinence, chronic yeast intertrigo, intestinal stomas, or in patients who desire to take deep tub baths.

Contraindications: This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy. An upper chest (presteral) exit-site location of the X-Series Catheter may be contraindicated in patients who have had breast augmentation or breast reconstruction or have a tracheostomy. However, this group of patients may be suitable candidates for an upper abdominal exit-site location of the X-Series Catheter.

BARBED LUER LOCK ADAPTOR

Indications for Use: The Medcomp® Barbed Luer Lock Adaptor is indicated for long-term use as a connection point for peritoneal dialysis catheters and male luers.

Contraindications: Use only with catheter provided by Medcomp®. Known or suspected allergies to the device.

TITANIUM BARBED LUER LOCK ADAPTOR

Indications for Use: The Medcomp® Titanium Barbed Luer Lock Adaptor is indicated for long-term use as a connection point for peritoneal dialysis catheters and male luers.

Contraindications: Use only with catheter provided by Medcomp®. Known or suspected allergies to the device.

PD CATHETER STYLET

Indications for Use: The Medcomp® Peritoneal Catheter Stylet is indicated for transient use to keep a peritoneal dialysis catheter straight and maneuverable during insertion into the peritoneum.

Contraindications: Use only with catheter provided by Medcomp®.

PD CATHETER TUNNELER

Indications for Use: The Medcomp® Peritoneal Catheter Tunneler is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.

X-SERIES CURVED TUNNELER

Indications for Use: The Medcomp® X-Series Peritoneal Catheter Curved Tunneler is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.

X-SERIES TUNNELING ROD


Indications for Use: The Medcomp® X-Series Peritoneal Catheter Tunneling Rod is indicated for transient use to dissect tissue to create the subcutaneous tunnel opening.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

V-Series



I-Series



X-Series



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcomp.net

ISO13485
Certified Quality System

PN2744GL REV. A 4/24



www.medcomp.net

I-Series

Material
Silicone

SIL

Tip Design
Straight & Coiled



French Size
15F



Cuffs
Single & Dual



Code	Description	Box Qty	Distal Tip to Cuff (cm)	Cuff to Cuff (cm)	Cuff to Coil (cm)
CATHETER ONLY STRAIGHT					
MPP41SGL	15F X 41CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	15	-	-
MPP46SGL	15F X 46CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	20	-	-
MSC31SGL	15F X 31CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	10.25	-	-
MSC37SGL	15F X 37CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	16.25	-	-
MSC42SGL	15F X 42CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	21	-	-
MSC46SGL	15F X 46CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	25	-	-
MDC31SGL	15F X 31CM I-SERIES CATHETER, DUAL CUFF	5	6.5	3	-
MDC32SGL	15F X 32.25CM I-SERIES CATHETER, DUAL CUFF	5	7	3.5	-
MDC37SGL	15F X 37CM I-SERIES CATHETER, DUAL CUFF	5	10	5	-
MDC42SGL	15F X 42CM I-SERIES CATHETER, DUAL CUFF	5	15	5	-
MDC47SGL	15F X 47CM I-SERIES CATHETER, DUAL CUFF	5	20	5	-
CATHETER ONLY COILED					
MPP40CGL	15F X 40.25CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	16.36	-	13
MPP57CGL	15F X 57.5CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	13	-	7
MPP60CGL	15F X 60CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	15	-	10
MDC57CGL	15F X 57.5CM I-SERIES CATHETER, DUAL CUFF	5	12	6	8
MDC63CGL	15F X 63CM I-SERIES CATHETER, DUAL CUFF	5	19	6	14
BASIC SET STRAIGHT					
MPD-146GL	15F X 46CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	20	-	-
MPD-146SGL	15F X 46CM I-SERIES CATHETER, SINGLE SUBCUTANEOUS CUFF	5	25	-	-
MPD-242GL	15F X 42CM I-SERIES CATHETER, DUAL CUFF	5	15	5	-
BASIC SET COILED					
MPD-160GL	15F X 60CM I-SERIES CATHETER, SINGLE PREPERITONEAL CUFF	5	34	-	15
MPD-257GL	15F X 57.5CM I-SERIES CATHETER, DUAL CUFF	5	31	6	12
MPD-263GL	15F X 63CM I-SERIES CATHETER, DUAL CUFF	5	38.25	6	19

Catheter Only Contents: (1) Catheter (1) Female Luer (1) Clamp (1) End Cap

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (6) Gauze (1) End Cap

V-Series

Material
Silicone

SIL

Tip Design
Straight & Coiled



French Size
15F



Cuffs
Dual



Code	Description	Box Qty	Distal Tip to Cuff (cm)	Cuff to Cuff (cm)	Cuff to Coil (cm)
CATHETER ONLY STRAIGHT					
MC200301GL	15F X 43CM V-SERIES CATHETER, LEFT	5	15	5	-
MC200302GL	15F X 43CM V-SERIES CATHETER, RIGHT	5	15	5	-
CATHETER ONLY COILED					
MC200401GL	15F X 63CM V-SERIES CATHETER, LEFT	5	17.82	5	13
MC200402GL	15F X 63CM V-SERIES CATHETER, RIGHT	5	17.82	5	13
MC200501GL	15F X 57CM V-SERIES CATHETER, LEFT	5	14.10	5	9
BASIC SET STRAIGHT					
MC200101GL	15F X 43CM V-SERIES CATHETER, LEFT	5	15	5	-
MC200102GL	15F X 43CM V-SERIES CATHETER, RIGHT	5	15	5	-
BASIC SET COILED					
MC200201GL	15F X 63CM V-SERIES CATHETER, LEFT	5	17.82	5	13
MC200202GL	15F X 63CM V-SERIES CATHETER, RIGHT	5	17.82	5	13
MC200502GL	15F X 57CM V-SERIES CATHETER, LEFT	5	14.10	5	9

Catheter Only Contents: (1) Catheter (1) Female Luer (1) Clamp (1) End Cap

Basic Set Contents: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Peelable Introducer (1) Female Luer (1) Clamp (6) Gauze (1) End Cap

X-Series

Material
Silicone

SIL

Tip Design
Coiled



French Size
15F



Code	Description	Box Qty	Cuff to Coil (cm)–Lower Catheter	Proximal Cuff to End of Catheter (cm)–Upper Catheter
CATHETER ONLY				
MC20XCGL	15F X 64CM X-SERIES CATHETER	5	12.75	28.5
BASIC SET				
MC20XSGL	15F X 64CM X-SERIES CATHETER	5	12.75	28.5

Catheter Only Contents: (1) Catheter (1) Titanium Barbed Connector (1) Female Luer (1) Clamp (1) End Cap

Basic Set Contents: (1) Catheter (1) Titanium Barbed Connector (1) Introducer Needle (1) Guidewire (1) Guidewire Advancer (1) Syringe (1) Scalpel (1) Tunneler (1) Female Luer (1) Clamp (6) Gauze (1) End Cap (1) Tape Measure