

Vascu-PICC TLS

PERIPHERALLY INSERTED CENTRAL CATHETER (PICC)

FEATURES & BENEFITS

LENGTHS

20cm, 50cm

LUMENS

Single & Dual


AVAILABLE IN

1.9F & 2.6F

Indicated for the administration of fluids, medications and nutritional therapy.

VERSATILITY

Since many pediatric patients may require a combination of nutrition, fluids and/or medicine, the dual lumen configuration enables one catheter to serve multiple functions.



Dual Lumen Design

Combined access to one site and offers the ability to infuse incompatible drugs

Vascu-PICC® TLS

PERIPHERALLY INSERTED CENTRAL CATHETER (PICC)

CATHETER ONLY SETS

MD17011100	1.9F X 50CM SINGLE LUMEN VASCU-PICC® TLS, CATHETER ONLY	10/BOX
MD17012600	2.6F X 50CM DUAL LUMEN VASCU-PICC® TLS, CATHETER ONLY	10/BOX

Catheter Only Set Contents: (1) PICC w/Stylet (1) Sideport Adaptor (1I2) Needleless Connector (1) Catheter Securement Device

BASIC SETS WITH OTN

MD17011121	1.9F X 20CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITH OTN	5/BOX
MD17011101	1.9F X 50CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITH OTN	5/BOX
MD17012621	2.6F X 20CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITH OTN	5/BOX
MD17012601	2.6F X 50CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITH OTN	5/BOX

Basic Set With OTN Contents: (1) PICC w/Stylet (1) Sideport Adaptor (1) Latex Free Tourniquet (3) 2" x 2" Gauze (1I2) Needleless Connector(s) (1) 10CC Syringe (1) Tape Measure (1) OTN Introducer (1) Catheter Securement Device

BASIC SETS WITH MST

MD170111224S	1.9F X 20CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITH MST	5/BOX
MD170111024S	1.9F X 50CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITH MST	5/BOX
MD170126224S	2.6F X 20CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITH MST	5/BOX
MD170126024S	2.6F X 50CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITH MST	5/BOX

Basic Set With MST Contents: (1) PICC w/Stylet (1) Sideport Adaptor (1) Latex Free Tourniquet (3) 2"x 2" Gauze (1I2) Needleless Connector(s) (1) 10CC Syringe (1) Tape Measure (1) Peelable Introducer (1) Guidewire (1) I.V. Catheter (1) Catheter Securement Device

BASIC SETS WITHOUT OTN

MD17021121	1.9F X 20CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITHOUT OTN	5/BOX
MD17021101	1.9F X 50CM SINGLE LUMEN VASCU-PICC® TLS, BASIC SET WITHOUT OTN	5/BOX
MD17022621	2.6F X 20CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITHOUT OTN	5/BOX
MD17022601	2.6F X 50CM DUAL LUMEN VASCU-PICC® TLS, BASIC SET WITHOUT OTN	5/BOX

Basic Set Without OTN Contents: (1) PICC w/Stylet (1) Sideport Adaptor (1) Latex Free Tourniquet (3) 2"x 2" Gauze (1I2) Needleless Connector(s) (1) 10CC Syringe (1) Tape Measure (1) Scalpel* (1) Catheter Securement Device

*2.6F Sets Only

FULL NURSING KIT W/MAX BARRIER CONTENTS

MRVSP19001	1.9F X 20CM SINGLE LUMEN VASCU-PICC® TLS, FULL NURSING KIT W/MAX BARRIER	3/BOX
MRVSP26001	2.6F X 50CM DUAL LUMEN VASCU-PICC® TLS, FULL NURSING KIT W/MAX BARRIER	3/BOX

Full Nursing Kit W/Max Barrier Contents: (1) PICC W/Stylet (1) Sideport Adaptor (1) 22ga Introducer Needle (1) I.V. Catheter (1) Guidewire (1) Scalpel (1) Peelable Introducer (1) 5CC Syringe (1) 25ga Safety Needle (1) Latex Free Tourniquet (3) Pre-Filled Saline Syringes (3) 2" x 2" Gauze (3) 4" x 4" Gauze (1) Lidocaine Ampule (1) Filter Straw (1) Probe Cover W/Gel (1) Towel (1) Gown (1) Fenestrated Drape (1) Drape (3) Clear Chlorhexidine (1) Scissors (1) Forceps (1) Alcohol Swabs (2) Face Masks (2) Tape Measures (1) Steri-Strips (1) Bouffant Cap (1) Catheter Securement Device (1) Needle Stick Pad (1) Needle Securement Device

PROCEDURE TRAY

MR17019000	PROCEDURE TRAY	3/BOX
------------	----------------	-------

Procedure Tray Set Contents: (3) Pre-Filled Saline Syringes (2) Chlorhexidine (1) Gown (1) Lidocaine Ampule (1) Gloves (2) Drapes (1) Fenestrated Drape (2) Face Masks (1) 5CC Syringe (1) 25ga Safety Needle (1) Scalpel (1) Adhesive Dressing (1) Filter Straw (1) Steri-Strips (1) Bouffant Cap (1) Needle Stick Pad (1) Scissors (2) Tape Measures (1) Lidocaine/Saline Stickers (3) 2" x 2" Gauze (3) 4" x 4" Gauze (1) Surgical Tape

PEDIATRIC PROCEDURE TRAY

MR17019001	PEDIATRIC PROCEDURE TRAY	3/BOX
------------	--------------------------	-------

Pediatric Procedure Tray Set Contents: (3) Pre-Filled Saline Syringes (2) Chlorhexidine (1) Gown (1) Lidocaine Ampule (1) Gloves (1) Drape (1) Fenestrated Drape (2) Face Masks (1) Steel Forceps (1) 5CC Syringe (1) 25ga Safety Needle (1) Scalpel (1) Adhesive Dressing (1) Filter Straw (1) Steri-Strips (1) Bouffant Cap (1) Needle Stick Pad (1) Scissors (2) Tape Measures (1) Lidocaine/Saline Stickers (3) 2" x 2" Gauze (3) 4" x 4" Gauze (1) Surgical Tape

IMPORTANT RISK INFORMATION

Indications for Use: The 1.9F Vascu-PICC® is a long-term catheter intended for central venous access via peripheral insertion in neonates, infants and children. It may be used for administration of fluids, medication and nutritional therapy.

Indications for Use: The 2.6F Vascu-PICC® catheters are indicated for short or long term access to the central venous system via peripheral insertion in neonates, infants and children. It may be used for the administration of fluids, medication and nutritional therapy.

Contraindications:

This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels
The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.)
The presence of device related bacteremia or septicemia.
Previous history of venous/subclavian thrombosis or vascular surgical procedures at insertion site.
Fever of unknown origin.
The patient's body size is insufficient to accommodate the size of the implanted device.
The patient is known or is suspected to be allergic to materials contained in the device.
Past irradiation of prospective insertion site.
Local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.



Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201
www.medcomp.net

ISO13485
Certified Quality System

PN2782US Rev. A 11/25