

# PD CATH

Peritoneal Dialysis Catheter

## FEATURES & BENEFITS

### TIP DESIGN

Straight & Coiled

### MATERIAL


Silicone

### AVAILABLE IN

15F

### CUFFS


Single & Double Cuff



V-Series



I-Series



Barbed Connector



X-Series



# PD CATH

## Peritoneal Dialysis Catheter

"X" SERIES		DESCRIPTION
<b>CATHETER ONLY 5/BOX</b>	<b>SETS 5/BOX</b>	
DMC20XC	DMC20XS	15F X 64CM X-SERIES
"V" SERIES		DESCRIPTION
<b>CATHETER ONLY 5/BOX</b>	<b>SETS 5/BOX</b>	
MC20VS43LC	MC20VS43LS	STRAIGHT, 43CM LEFT V-NECK PD CATHETER WITH DOUBLE CUFF
MC20VS43RC	MC20VS43RS	STRAIGHT, 43CM RIGHT V-NECK PD CATHETER WITH DOUBLE CUFF
MC20VC57LC	MC20VC57LS	COILED, 57CM LEFT V-NECK PD CATHETER WITH DOUBLE CUFF
MC20VC63LC	MC20VC63LS	COILED, 63CM LEFT V-NECK PD CATHETER WITH DOUBLE CUFF
MC20VC63RC	MC20VC63RS	COILED, 63CM RIGHT V-NECK PD CATHETER WITH DOUBLE CUFF
"I" SERIES		DESCRIPTION
<b>CATHETER ONLY 5/BOX</b>	<b>SETS 5/BOX</b>	
MC20IS31DC	—	STRAIGHT, 31CM PD CATHETER WITH DOUBLE CUFF
MC20IS32DC	—	STRAIGHT, 32.25CM PD CATHETER WITH DOUBLE CUFF
MC20IS37DC	—	STRAIGHT, 37CM PD CATHETER WITH DOUBLE CUFF
MC20IS42DC	MC20IS42DS	STRAIGHT, 42CM PD CATHETER WITH DOUBLE CUFF
MC20IS47DC	—	STRAIGHT, 47CM PD CATHETER WITH DOUBLE CUFF
MC20IC57DC	MC20IC57DS	COILED, 57.5 PD CATHETER WITH DOUBLE CUFF
MC20IC63DC	MC20IC63DS	COILED, 63CM PD CATHETER WITH DOUBLE CUFF
MC20IS31SC	—	STRAIGHT, 31CM PD CATHETER WITH SUBCUTANEOUS SINGLE CUFF
MC20IS37SC	—	STRAIGHT, 37CM PD CATHETER WITH SUBCUTANEOUS SINGLE CUFF
MC20IS42SC	—	STRAIGHT, 42CM PD CATHETER WITH SUBCUTANEOUS SINGLE CUFF
MC20IS46SC	MC20IS46SS	STRAIGHT, 46CM PD CATHETER WITH SUBCUTANEOUS SINGLE CUFF
MC20IS41PC	—	STRAIGHT, 41CM PD CATHETER WITH PREPERITONEAL SINGLE CUFF
MC20IS46PC	MC20IS46PS	STRAIGHT, 46CM PD CATHETER WITH PREPERITONEAL SINGLE CUFF
MC20IC40PC	—	COILED, 40.25CM PD CATHETER WITH PREPERITONEAL SINGLE CUFF
MC20IC57PC	—	COILED, 57.5CM PD CATHETER WITH PREPERITONEAL SINGLE CUFF
MC20IC60PC	MC20IC60PS	COILED, 60CM PD CATHETER WITH PREPERITONEAL SINGLE CUFF

**Catheter Only Set Components:** (1) Catheter (1) Female Luer (1) Clamp (1) Male Luer End Cap

**Set Components:** (1) Catheter (1) Female Luer (1) Clamp (1) Male Luer End Cap (1) Introducer Dilator & Sheath (1) J/Flex Guidewire (1) Scalpel (1) Introducer Needle (1) Tunneling Tool (6) Gauze Pads (1) 10cc Luer Lock Syringe

ACCESSORIES		QTY
DMP217	17F PEELABLE INTRODUCER SET	5/BOX
DMBCA	BARBED LUER LOCK CATHETER ADAPTOR	20/BOX
DMCS	CATHETER STYLET	5/BOX
DMTS	TUNNELING STYLET	10/BOX
DMTC	TITANIUM CONNECTOR LUER	5/BOX
DTBC	TITANIUM BARBED CONNECTOR	50/BOX
DFSX	X-SERIES FALLER STYLET	5/BOX
DTRX	X-SERIES TUNNELING ROD	5/BOX
40501	PREOPERATIVE MARKING STENCIL PACK	1/EACH

### IMPORTANT RISK INFORMATION

#### Peritoneal Dialysis (I-Series, V-Series, X-Series)

**Indications for Use:** The Peritoneal Catheter is indicated for chronic peritoneal dialysis. The extended length of the "X" Series Catheter makes it especially applicable for peritoneal dialysis patients when it is necessary to locate the skin exit site remote from the usual lower abdominal region. The catheter may be particularly indicated in patients with obesity, floppy abdominal skin folds, urinary or fecal incontinence, chronic yeast intertrigo, intestinal stomas, or in patients who desire to take deep tub baths

**Contraindications:** This catheter is contraindicated in patients who are not suitable candidates for peritoneal dialysis therapy. An upper chest (presternal) exit-site location of the X-Series Catheter may be contraindicated in patients who have had breast augmentation or breast reconstruction, or have a tracheostomy. However, this group of patients may be suitable candidates for an upper abdominal exit-site location of the X-Series Catheter.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.



**Medical Components, Inc.**  
1499 Delp Drive, Harleysville, PA 19438, U.S.A.  
P: 215.256.4201 F: 215.256.1787  
www.medcompnet.com

**ISO13485**  
Certified Quality System

PN2453 REV. I 6/21