

14F & 16F CATHETER SET PRE-CURVED

ASPC24-3PCE.	14F X 24CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC28-3PCE.	14F X 28CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC32-3PCE.	14F X 32CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC36-3PCE.	14F X 36CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC2416-3PCE.	16F X 24CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC2816-3PCE.	16F X 28CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC3216-3PCE.	16F X 32CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX
ASPC3616-3PCE.	16F X 36CM SPLIT CATH® III CATHETER, PRE-CURVED	5/BOX

14F & 16F PRE-CURVED CATHETER SET CONTENTS: (1) Catheter (1) Introducer Needle (1) Guidewire with Advancer** (1) Dilator (1) Tunneler (1) Valved Peelable Introducer (2) End Cap (1) Patient Information Packet (1) Patient ID Card

** — 32 & 36CM lengths include (2) Guidewires With Advancer

14F CATHETER SET WITHOUT SIDEHOLE

ASPC24-3WOE.	14F X 24CM SPLIT CATH® III CATHETER WITHOUT SIDEHOLE, STRAIGHT	5/BOX
ASPC28-3WOE.	14F X 28CM SPLIT CATH® III CATHETER WITHOUT SIDEHOLE, STRAIGHT	5/BOX
ASPC32-3WOE.	14F X 32CM SPLIT CATH® III CATHETER WITHOUT SIDEHOLE, STRAIGHT	5/BOX
ASPC36-3WOE.	14F X 36CM SPLIT CATH® III CATHETER WITHOUT SIDEHOLE, STRAIGHT	5/BOX
ASPC40-3WOE.	14F X 40CM SPLIT CATH® III CATHETER WITHOUT SIDEHOLE, STRAIGHT	5/BOX

14F WITHOUT SIDEHOLE CATHETER SET CONTENTS: (1) Catheter With Stylet (1) Introducer Needle (1) Guidewire with Advancer** (1) Dilator (1) Tunneler (1) Valved Peelable Introducer (2) End Cap (1) Patient Information Packet (1) Patient ID Card

** — 32, 36 & 40CM lengths include (2) Guidewires With Advancer

14F TRANSLUMBAR SET

ASPC40-3TLE.	14F X 40CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC55-3TLE.	14F X 55CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX

14F TRANSLUMBAR SET CONTENTS: (1) Catheter With Stylet (1) Introducer Needle (2) Guidewire with Advancer (3) Dilator (1) Tunneler (1) Peelable Introducer (1) Scalpel (2) End Cap (1) Patient Information Packet (1) Patient ID Card

IMPORTANT RISK INFORMATION

Split Cath®

Indications for Use: The Split Cath® Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis.

Split Cath® III

Indications for Use: The Split Cath® III Catheter is indicated for short-term or long-term use where vascular access is required for 14 days or more for the purpose of hemodialysis and apheresis.

Split Cath® & Split Cath® III

Contraindications: Known or suspected allergies to any of the components of the catheter or the kit. This device is contraindicated for patients exhibiting severe, uncontrolled coagulopathy or thrombocytopenia.

Refer to Instructions for Use provided with the product for complete instructions, warnings, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

SPLIT CATH® & SPLIT CATH® III

LONG-TERM HEMODIALYSIS CATHETER



FEATURES & BENEFITS

LUMEN DESIGN

Double “D” Internal



TIP DESIGN

Split Tip



AVAILABLE IN

10F, 14F & 16F



CONFIGURATIONS

Straight



Pre-Curved



SPLIT CATH®

LONG-TERM HEMODIALYSIS CATHETER

10F CATHETER SET

ASPC15P-XL	10F X 15CM SPLIT CATH® CATHATER, STRAIGHT	5/BOX
ASPC18P-XL	10F X 18CM SPLIT CATH® CATHATER, STRAIGHT	5/BOX
ASPC24P-XL	10F X 24CM SPLIT CATH® CATHATER, STRAIGHT	5/BOX

10F CATHETER SET CONTENTS: (1) Catheter (1) Introducer Needle (1) Guidewire With Advancer (2) Dilator (1) Tunneler (1) Peelable Introducer (1) Scalpel (2) End Cap (1) Patient Information Packet (1) Patient ID Card

SPLIT CATH® III

LONG-TERM HEMODIALYSIS CATHETER

14F & 16F STRAIGHT CATHETER SET

ASPC20-3E.	14F X 20CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC22-3E.	14F X 22CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC24-3	14F X 24CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC24-3E.	14F X 24CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC28-3	14F X 28CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC28-3E.	14F X 28CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC32-3	14F X 32CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC32-3E.	14F X 32CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC36-3	14F X 36CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC36-3E.	14F X 36CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC40-3	14F X 40CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC40-3E.	14F X 40CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC55-3	14F X 55CM SPLIT CATH® III CATHETER, STRAIGHT*	5/BOX
ASPC55-3E.	14F X 55CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC2416-3E.	16F X 24CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC2816-3E.	16F X 28CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC3216-3E.	16F X 32CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC3616-3E.	16F X 36CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC4016-3E.	16F X 40CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX
ASPC5516-3E.	16F X 55CM SPLIT CATH® III CATHETER, STRAIGHT	5/BOX

14F & 16F STRAIGHT CATHETER SET CONTENTS: (1) Catheter With Stylet (1) Introducer Needle (1) Guidewire With Advancer** (1) Dilator (1) Tunneler (1) Valved Peelable Introducer (2) End Cap (1) Patient Information Packet (1) Patient ID Card

* — Codes include Scalpel ** — 32, 36 & 40CM lengths include (2) Guidewires With Advancer