



# medCOMP<sup>®</sup> MST Kit



## Modified Seldinger Technique



### IMPORTANT RISK INFORMATION

#### MODIFIED SELDINGER TECHNIQUE (MST) KIT


**Indications for Use:** The 2F and 3F Vascu-Sheath Tearaway Introducer is intended for percutaneous venous access by modified Seldinger Technique in neonates, infants, and children.

#### Contraindications:

- This catheter is not intended for any use other than that which is indicated. Do not implant catheter in thrombosed vessels.
- The presence of skin related problems around the insertion site (infection, phlebitis, scars, etc.)
- The presence of device related bacteremia or septicemia.
- Previous history of venous/ subclavian thrombosis or vascular surgical procedures at insertion site.
- Fever of unknown origin.
- The patient's body size is insufficient to accommodate the size of the implanted device.
- The patient is known or is suspected to be allergic to materials contained in the device.
- Past irradiation of prospective insertion site.
- Local tissue factors will prevent proper device stabilization and/ or access.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.

# medCOMP<sup>®</sup> MST Kit



2F MST KITS		UOM
DMR1902224S	2F MST W/ 20CM WIRE & 24GA SAFETY	10/BOX
DMR1902222S	2F MST W/ 20CM WIRE & 22GA SAFETY	10/BOX
DMR1902222	2F MST W/ 20CM WIRE & 22GA	10/BOX

2F MST KITS CONTAINS:  
 (1) Tearaway Introducer  
 (1) Needle  
 (1) Guidewire

DMR1902424S	2F MST W/ 40CM WIRE & 24GA SAFETY	10/BOX
DMR1902422S	2F MST W/ 40CM WIRE & 22GA SAFETY	10/BOX
DMR1902422	2F MST W/ 40CM WIRE & 22GA	10/BOX

3F MST KITS CONTAINS:  
 (1) Tearaway Introducer  
 (1) Needle  
 (1) Guidewire

DMR1902624S	2F MST W/ 60CM WIRE & 24GA SAFETY	10/BOX
DMR1902622S	2F MST W/ 60CM WIRE & 22GA SAFETY	10/BOX
DMR1902622	2F MST W/ 60CM WIRE & 22GA	10/BOX

3F MST KITS		UOM
DMR1903224S	3F MST W/ 20CM WIRE & 24GA SAFETY	10/BOX
DMR1903222S	3F MST W/ 20CM WIRE & 22GA SAFETY	10/BOX
DMR1903222	3F MST W/ 20CM WIRE & 22GA	10/BOX

DMR1903424S	3F MST W/ 40CM WIRE & 24GA SAFETY	10/BOX
DMR1903422S	3F MST W/ 40CM WIRE & 22GA SAFETY	10/BOX
DMR1903422	3F MST W/ 40CM WIRE & 22GA	10/BOX

DMR1903624S	3F MST W/ 60CM WIRE & 24GA SAFETY	10/BOX
DMR1903622S	3F MST W/ 60CM WIRE & 22GA SAFETY	10/BOX
DMR1903622	3F MST W/ 60CM WIRE & 22GA	10/BOX

ACCESSORIES		UOM
DVS203	2F TEARAWAY INTRODUCER	5/BOX
DVS303	3F TEARAWAY INTRODUCER	5/BOX
MR190020	.010 X 20CM NITINOL GUIDE WIRE	10/BOX
DMR190040	.010 X 40CM NITINOL GUIDE WIRE	10/BOX
DMR190060	.010 X 60CM NITINOL GUIDE WIRE	10/BOX
DMR19100801	.010 X 80CM NITINOL GUIDE WIRE	10/BOX