



Important Risk Information

Indications For Use: The Pro-Fuse® Power Injectable Port is indicated for long-term access to the central venous system for intravenous administration of fluids or medications, power injection of contrast media, and withdrawal of blood samples.

Contraindications: This device is contraindicated for catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates of pinch-off.

The device is also contraindicated:

- When the presence of device related infection, bacteremia, or septicemia is known or suspected.
- When the patient's body size is insufficient for the size of the implanted device.
- When the patient is known or is suspected to be allergic to materials contained in the device.
- If severe chronic obstructive lung disease exists.
- If the prospective insertion site has been previously irradiated.
- If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures.
- If local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation



SIZES

Pro-Fuse®
Low Profile



Pro-Fuse®



FEATURES & BENEFITS



Universal Locking Collar



CT Printed in reverse
on port base.



Radiopaque Ink



medcomp.net

ORDERING INFORMATION

	OPEN SUTURE HOLES	SILICONE FILLED SUTURE HOLES	
STANDARD SET			
8F Attachable	–	MRCTT80001	5/BOX
9.6F Attachable	–	MRCTT96001	5/BOX
8F Pre-Attached	–	MRCTT80801	5/BOX
9.6F Pre-Attached**	–	MRCTT96801	5/BOX
LOW PROFILE SET			
6.6F Attachable	MRCTT66041	MRCTT6604S	5/BOX
8F Attachable	MRCTT80041	–	5/BOX
6.6F Pre-Attached	MRCTT66841	–	5/BOX
STANDARD DMP SET			
8F Attachable	–	MRCTT80001DMP	5/BOX
LOW PROFILE DMP SET			
8F Attachable	–	MRCTT800041DMP	5/BOX

PRO-FUSE® POWER INJECTABLE PORT SET CONTENTS: (1) Port (1) Catheter (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Guidewire (1) Advancer (1) Syringe (1) Peelable Introducer* (1) Tunneler (2) Huber Needles (1) Blunt Tip Needle (1) Vein Pick (1) Patient Information Pack (1) Patient ID Card

PRO-FUSE® POWER INJECTABLE PORT WITH DIRECT MICROPUNCTURE SET CONTENTS: (1) Port (1) Catheter (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Guidewire (1) Syringe (1) Peelable Introducer* (1) Tunneler (1) Huber Needle (1) Blunt Tip Needle (1) Patient Information Pack (1) Patient ID Card

*6.6F Kit has a 7F Peelable Introducer

8F Kit has a 9F Peelable Introducer

9.6F Kit has a 10F Peelable Introducer

**9.6F Kit has a Silicone Catheter