



Important Risk Information

Indications For Use: The Dignity[®] Power Injectable Port is indicated for long-term access to the central venous system for intravenous administration of fluids or medications, power injection of contrast media, and withdrawal of blood samples.

Contraindications: This device is contraindicated for catheter insertion in the subclavian vein medial to the border of the first rib, an area which is associated with higher rates of pinch-off.

The device is also contraindicated:

- When the presence of device related infection, bacteremia, or septicemia is known or suspected.
- When the patient's body size is insufficient for the size of the implanted device.
- When the patient is known or is suspected to be allergic to materials contained in the device.
- If severe chronic obstructive lung disease exists.
- If the prospective insertion site has been previously irradiated.
- If the prospective placement site has previously suffered episodes of venous thrombosis or vascular surgical procedures.
- If local tissue factors will prevent proper device stabilization and/or access.

Refer to Instructions for Use for complete instructions, warnings, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.

For patient information please visit www.medcomp.net/patientinformation

Dignity®

Power Injectable Ports

SIZES

Low Profile



Mid Sized



Mini



FEATURES & BENEFITS

Tapered Tail



Radiopaque Ink



Comfort



Reduced Width



Universal Locking Collar



Smooth Contours



Raised Rim

Sculpted Sides

ORDERING INFORMATION

	OPEN SUTURE HOLES		SILICONE FILLED SUTURE HOLES		
		WITH MICRO-STICK®		WITH MICRO-STICK®	
MINI PROFILE SET					
5F ATTACHABLE	MRCTI50041M	MICTI50041M	MRCTI5004SM	MICTI5004SM	5/BOX
6.6F ATTACHABLE	MRCTI66041M	MICTI66041M	MRCTI6604SM	MICTI6604SM	5/BOX
8F ATTACHABLE	MRCTI80041M	MICTI80041M	-	MICTI8004SM	5/BOX
5F PRE-ATTACHED	MRCTI50841M	-	MRCTI5084SM	-	5/BOX
6.6F PRE-ATTACHED	-	-	MRCTI6684SM	-	5/BOX
8F PRE-ATTACHED	-	-	-	MICTI8084SM	5/BOX

LOW PROFILE SET					
5F ATTACHABLE	MRCTI50041	-	-	MICTI5004S	5/BOX
6.6F ATTACHABLE	MRCTI66041	MICTI66041	MRCTI6604S	MICTI6604S	5/BOX
8F ATTACHABLE	MRCTI80041	MICTI80041	MRCTI8004S	MICTI8004S	5/BOX
6.6F PRE-ATTACHED	MRCTI66841	MICTI66841	MRCTI6684S	-	5/BOX
8F PRE-ATTACHED	MRCTI80841	MICTI80841	MRCTI8084S	-	5/BOX

MID-SIZED SET					
5F ATTACHABLE	MRCTI50001	-	-	-	5/BOX
6.6F ATTACHABLE	MRCTI66001	MICTI66001	MRCTI6600S	MICTI6600S	5/BOX
8F ATTACHABLE	MRCTI80001	MICTI80001	MRCTI8000S	MICTI8000S	5/BOX
9.6F ATTACHABLE**	MRCTI96001	-	MRCTI9600S	-	5/BOX
6.6F PRE-ATTACHED	MRCTI66801	-	-	-	5/BOX
8F PRE-ATTACHED	MRCTI80801	-	MRCTI8080S	-	5/BOX
9.6F PRE-ATTACHED**	MRCTI96801	-	-	-	5/BOX

MINI PROFILE DMP SET					
5F ATTACHABLE DMP	MRCTI50041MDMP	-	-	-	5/BOX
6.6F ATTACHABLE DMP	MRCTI66041MDMP	-	-	-	5/BOX
8F ATTACHABLE DMP	MRCTI80041MDMP	-	-	-	5/BOX

LOW PROFILE DMP SET					
5F ATTACHABLE DMP	MRCTI50041DMP	-	-	-	5/BOX
6.6F ATTACHABLE DMP	MRCTI66041DMP	-	-	-	5/BOX
8F ATTACHABLE DMP	MRCTI80041DMP	-	-	-	5/BOX

MID-SIZED DMP SET					
6.6F ATTACHABLE DMP	MRCTI66001DMP	-	-	-	5/BOX
8F ATTACHABLE DMP	MRCTI80001DMP	-	-	-	5/BOX
6.6F PRE-ATTACHED DMP	MRCTI66801DMP	-	-	-	5/BOX

DIGNITY® POWER INJECTABLE PORT SET CONTENTS: (1) Port (1) Catheter (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Guidewire (1) Advancer (1) Syringe (1) Peelable Introducer* (1) Tunneler (2) Huber Needles (1) Blunt Tip Needle (1) Vein Pick (1) Patient Information Pack (1) Patient ID Card

DIGNITY® POWER INJECTABLE PORT WITH MICRO-STICK® SET CONTENTS: (1) Port (1) Catheter (2) Catheter Locks (1) Scalpel (2) Guidewires (1) Advancer (1) Syringe (1) Peelable Introducer* (1) Tunneler (2) Huber Needles (1) Blunt Tip Needle (1) Vein Pick (1) Patient Information Pack (1) Patient ID Card (1) Micro-Stick® (1) Introducer Needle

DIGNITY® POWER INJECTABLE PORT WITH DIRECT MICROPUNCTURE SET CONTENTS: (1) Port (1) Catheter (2) Catheter Locks (1) Scalpel (1) Introducer Needle (1) Guidewire (1) Syringe (1) Peelable Introducer* (1) Tunneler (1) Huber Needle (1) Blunt Tip Needle (1) Patient Information Pack (1) Patient ID Card

*5F Kit has a 5.5F Peelable Introducer

6.6F Kit has a 7F Peelable Introducer

8F Kit has a 9F Peelable Introducer

9.6F Kit has a 10F Peelable Introducer

**9.6F Kit has a Silicone Catheter