



TESIO®

**Indications for Use:** The Medcomp® Tesio® Catheter is designed for long term hemodialysis and apheresis. It may be inserted percutaneously and is ideally placed in the internal jugular vein. Although this catheter may be inserted into the subclavian vein, the internal jugular is the preferred site.


**Contraindications:** This catheter is intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions. To maintain peak performance, it is recommended that this catheter be replaced after 18 months usage. To maintain peak performance of the Lock right Adapters, it is recommended that the adapters be replaced every 6 months.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.



<b>Polyurethane</b>	<b>Single Round</b>	<b>6.5F, 10F</b>	<b>Straight</b>	<b>Retrograde</b>
Material	Tip Design	French Sizes	Configuration	Insertion Type

  
**INDEPENDENT SINGLE  
 ROUND INTERNAL  
 LUMEN DESIGN**



**6.5F BASIC SET**

PBF	6.5F X 29CM TESIO® CATHETER, CUFFS @ 12CM AND 15CM	5/BOX
-----	--	-------

**10F BASIC SET**

DBFS-6	10F X 52CM TESIO® CATHETERS, CUFFS @ 18CM AND 21CM	5/BOX
DBFR-6	10F X 52CM TESIO® CATHETERS, CUFFS @ 22CM AND 25CM	5/BOX
DBFL-6	10F X 52CM TESIO® CATHETERS, CUFFS @ 27CM AND 30CM	5/BOX
BFRD1070KD-A	10F X 70CM TESIO® CATHETERS, CUFFS @ 46CM AND 50CM	5/BOX

**52CM BASIC SET**

DBFR1035KD-A	10F X 52CM TESIO® CATHETERS, CUFFS @ 18CM AND 22CM	5/BOX
--------------	--	-------

**BASIC SET WITH STYLET**

DBFS-6S*	10F X 52CM TESIO® CATHETERS, CUFFS @ 18CM AND 21CM	5/BOX
DBFR-6S*	10F X 52CM TESIO® CATHETERS, CUFFS @ 22CM AND 25CM	5/BOX
DBFL-6S*	10F X 52CM TESIO® CATHETERS, CUFFS @ 27CM AND 30CM	5/BOX

**TESIO® REPAIR KIT**

PBF11A	6.5F TESIO® CATHETER ARTERIAL EXTENSION SET	10/BOX
DBF11A	10F TESIO® CATHETER ARTERIAL EXTENSION SET	10/BOX
PBF11V	6.5F TESIO® CATHETER VENOUS EXTENSION SET	10/BOX
DBF11V	10F TESIO® CATHETER VENOUS EXTENSION SET	10/BOX
DBF11AVS	10F TESIO® VENOUS/ARTERIAL EXTENSION SETS	10/BOX

**6.5F BASIC SET CONTENTS:**

- (2) Catheters
- (2) Introducer Needles
- (2) Guidewires
- (2) Tunneler
- (2) Peelable Introducers
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (2) Clamps
- (2) End Caps
- (1) Patient ID Card

**10F BASIC SET CONTENTS:**

- (2) Catheters\*
- (2) Introducer Needles
- (2) Guidewires
- (2) Tunneler
- (2) Football Tunneler
- (1) Dilator
- (2) Valved Peelable Introducers
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (2) Clamps
- (2) End Caps
- (1) Securement Device
- (1) Patient ID Card

**52CM BASIC SET CONTENTS:**

- (2) Catheters
- (2) Introducer Needles
- (2) Guidewires
- (2) Tunneler
- (2) Valved Peelable Introducers
- (1) Arterial Extension Set (Red)
- (1) Venous Extension Set (Blue)
- (2) Clamps
- (2) End Caps
- (1) Patient ID Card

**TESIO® ARTERIAL EXTENSION SET CONTENTS:**

- (1) Arterial Extension Set
- (1) End Cap

**TESIO® VENOUS EXTENSION SET CONTENTS:**

- (1) Venous Extension Set
- (1) End Cap

**TESIO® VENOUS/ARTERIAL EXTENSION SET CONTENTS:**

- (1) Arterial Extension Set
- (1) Venous Extension Set
- (2) End Cap

\* Stylet Included

Refer to the Table of Contents for Important Risk Information regarding this device.