

ACCESSORIES

DRAINAGE INTRODUCER KITS

BINTN	BREEZE® INTRODUCER KIT WITH NITINOL & STAINLESS STEEL WIRE	5/BOX
BINTSS	BREEZE® INTRODUCER KIT WITH STAINLESS STEEL WIRE	5/BOX

BINTN Kit Contents: (1) 6F x 20CM Introducer with Hydrophillic Coating (1) 0.018" x 80CM Nitinol Guidewire with Radiopaque Tip (1) 0.035" x 150CM Stainless Steel Guidewire (1) 21GA x 15CM Trocar Introducer Needle, Echogenic

BINTSS Kit Contents: (1) 6F x 20CM Introducer with Hydrophillic Coating (1) 0.018" x 80CM Stainless Steel Guidewire with Radiopaque Tip (1) 0.035" x 150CM Stainless Steel Guidewire (1) 21GA x 15CM Trocar Introducer Needle, Echogenic



DRAINAGE SECUREMENT DEVICE

BL0512	BREEZE® DRAIN LOCK 5F-12F	25/BOX
BL1218	BREEZE® DRAIN LOCK 12F-18F	25/BOX



IMPORTANT RISK INFORMATION

BREEZE® and BREEZE® XS All Purpose Drain

Indications for Use: The catheter is designed for percutaneous drainage of abscess fluid, cyst, gall bladders, nephrostomy, urinary, and others fluids. The operator can use different drainage sets according to the type of accumulated fluid and place of accumulation. These sets are classified according to the catheter size (5F-16F pertaining to various fluid viscosity) and according to the position of the accumulation. The operator can BREEZE® and BREEZE® XS DRAINAGE CATHETER SET choose from either the Direct Access or Seldinger Type depending on the method used to provide access.

Contraindications: Only qualified personnel who are familiar with the technique should use the product. During insertion avoid contact with bone, cartilage and scar tissue which may damage the catheter tip. The set should not be used where percutaneous catheterization is unacceptable. Echinococcus cyst, multiple chamber accumulation of fluids or accumulations known as phlegmonic accumulations can mostly not be satisfactorily treated using the percutaneous drainage technique.

ALL PURPOSE DRAIN ACCESSORIES

Indications for Use: The Breeze® and Breeze® XS Introducer System is indicated for use in percutaneous introduction and placement of catheters and guidewires.

Contraindications: Use of the introducer is contraindicated if the patient has a known or suspected obstruction in the vessel or target anatomy. There is increased risk of pneumothorax for the patient who has severe chronic lung disease.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all Instructions for Use prior to using products. Failure to do so may result in patient complications.




Medical Components, Inc.
1499 Delp Drive, Harleysville, PA 19438, U.S.A.
P: 215.256.4201 F: 215.256.1787
www.medcompnet.com

ISO13485
Certified Quality System

PN2718 REV. A 10/21

BREEZE® / BREEZE® XS

ALL PURPOSE DRAINAGE CATHETER



medCOMP
www.medcompnet.com


BREEZE®

Drainage Catheter

Tip Design Pigtail with Sideholes

Tip Options Locking, Non-Locking & Locking Mini

Available In
5F, 6F, 7F, 8F, 10F, 12F, 14F & 16F



ALL PURPOSE DRAIN WITH LOCKING PIGTAIL			DIM A
BALPN0625S	6F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN0630S	6F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN0725S	7F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN0730S	7F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN0820S	8F X 20CM APD WITH LOCKING PIGTAIL	5/BOX	20CM
BALPN0825S	8F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN0830S	8F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN0835S	8F X 35CM APD WITH LOCKING PIGTAIL	5/BOX	35CM
BALPN1020S	10F X 20CM APD WITH LOCKING PIGTAIL	5/BOX	20CM
BALPN1025S	10F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN1030S	10F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN1035S	10F X 35CM APD WITH LOCKING PIGTAIL	5/BOX	35CM
BALPN1225S	12F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN1230S	12F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN1425S	14F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN1430S	14F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
BALPN1625S	16F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM
BALPN1630S	16F X 30CM APD WITH LOCKING PIGTAIL	5/BOX	30CM
ALL PURPOSE DRAIN WITHOUT LOCKING PIGTAIL			DIM A
BANPN0620S	6F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM
BANPN0625S	6F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN0720S	7F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM
BANPN0725S	7F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN0820S	8F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM
BANPN0825S	8F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN1020S	10F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM
BANPN1025S	10F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN1225S	12F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN1425S	14F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN1625S	16F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM
BANPN1625S	16F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM

ALL PUPOSE DRAIN WITH LOCKING MINI PIGTAIL			DIM A
BALMN0620S	6F X 20CM APD WITH LOCKING MINI PIGTAIL	5/BOX	20CM
BALMN0820S	8F X 20CM APD WITH LOCKING MINI PIGTAIL	5/BOX	20CM
BALMN0825S	8F X 25CM APD WITH LOCKING MINI PIGTAIL	5/BOX	25CM
BALMN1020S	10F X 20CM APD WITH LOCKING MINI PIGTAIL	5/BOX	20CM
BALMN1025S	10F X 25CM APD WITH LOCKING MINI PIGTAIL	5/BOX	25CM

6F Kit Contents: (1) Catheter (1) Metal Trocar Needle (1) Metal Stiffening Cannula

7F-16F Kit Contents: (1) Catheter (1) Metal Trocar Needle (1) Metal Stiffening Cannula (1) Flexible Stiffening Cannula


BREEZE® XS

Drainage Catheter

Tip Design Pigtail with Sideholes

Tip Options Locking, Non-Locking & Locking Mini

Available In
5F, 6F, 7F, 8F, 10F & 12F



ALL PURPOSE DRAIN WITH LOCKING PIGTAIL			DIM A	DIM B
XSALPN0620S	6F X 20CM APD WITH LOCKING PIGTAIL	5/BOX	20CM	10CM
XSALPN0625S	6F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM	15CM
XSALPN0720S	7F X 20CM APD WITH LOCKING PIGTAIL	5/BOX	20CM	10CM
XSALPN0725S	7F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM	15CM
XSALPN0825S	8F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM	15CM
XSALPN1020S	10F X 20CM APD WITH LOCKING PIGTAIL	5/BOX	20CM	10CM
XSALPN1025S	10F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM	15CM
XSALPN1225S	12F X 25CM APD WITH LOCKING PIGTAIL	5/BOX	25CM	15CM
ALL PURPOSE DRAIN WITHOUT LOCKING PIGTAIL			DIM A	DIM B
XSANPN0520S	5F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM	15CM
XSANPN0620S	6F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM	15CM
XSANPN0720S	7F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM	15CM
XSANPN0825S	8F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM	15CM
XSANPN1020S	10F X 20CM APD WITHOUT LOCKING PIGTAIL	5/BOX	20CM	10CM
XSANPN1025S	10F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM	15CM
XSANPN1225S	12F X 25CM APD WITHOUT LOCKING PIGTAIL	5/BOX	25CM	15CM
ALL PUPOSE DRAIN WITH LOCKING MINI PIGTAIL			DIM A	DIM B
XSALMN0620S	6F X 20CM APD WITH LOCKING MINI PIGTAIL	5/BOX	20CM	15CM
XSALMN0825S	8F X 25CM APD WITH LOCKING MINI PIGTAIL	5/BOX	25CM	15CM
XSALMN1025S	10F X 25CM APD WITH LOCKING MINI PIGTAIL	5/BOX	25CM	15CM

5F & 6F Kit Contents: (1) Catheter (1) Metal Trocar Needle (1) Metal Stiffening Cannula

7F-16F Kit Contents: (1) Catheter (1) Metal Trocar Needle (1) Metal Stiffening Cannula (1) Flexible Stiffening Cannula

CODES CONTINUED ➔