

TITAN^{HD}

LONG-TERM HEMODIALYSIS CATHETER



FEATURES & BENEFITS

INTERNAL LUMEN DESIGN

Double "D"



MATERIAL

Polyurethane

TIP DESIGN

Step Tip



FRENCH SIZES

15.5F



CONFIGURATIONS

Straight



Pre-Curved





WITH SIDHOLES

BASIC SET WITH STYLET

THD155024SCA	15.5F X 24CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155028SCA	15.5F X 28CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155032SCA	15.5F X 32CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155036SCA	15.5F X 36CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155040SCA	15.5F X 40CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155055SCA	15.5F X 55CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155224SCA	15.5F X 24CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155228SCA	15.5F X 28CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155232SCA	15.5F X 32CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155236SCA	15.5F X 36CM TITAN HD CATHETER, PRE-CURVED	5/BOX

Basic Set Contents w/ Sideholes: (1) Catheter (1) Stylet (1) Introducer Needle (1) Guidewire (1) Tunneler (2) Dilator (1) Valved Peelable Introducer (2) End Cap

BASIC SET WITHOUT STYLET

THD155224CA	15.5F X 24CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155228CA	15.5F X 28CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155232CA	15.5F X 32CM TITAN HD CATHETER, PRE-CURVED	5/BOX
THD155236CA	15.5F X 36CM TITAN HD CATHETER, PRE-CURVED	5/BOX

Basic Set Contents w/o Sideholes w/o Stylet: (1) Catheter (1) Introducer Needle (1) Guidewire (1) Tunneler (2) Dilator (1) Valved Peelable Introducer (2) End Cap

WITHOUT SIDHOLES

BASIC SET

THD155424SCA	15.5F X 24CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155428SCA	15.5F X 28CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155432SCA	15.5F X 32CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155436SCA	15.5F X 36CM TITAN HD CATHETER, STRAIGHT	5/BOX
THD155440SCA	15.5F X 40CM TITAN HD CATHETER, STRAIGHT	5/BOX

Basic Set Contents w/o Sideholes: (1) Catheter, (1) Stylet, (1) Introducer Needle (1) Guidewire (1) Tunneler (2) Dilator (1) Valved Peelable Introducer (2) End Cap

ACCESSORIES

RPK-01CA	MEDCOMP [®] REPAIR KIT	5/BOX
----------	---------------------------------	-------

Repair Kit Contents: (1) Luer Assembly (1) Scissors (1) Adhesive Backed Drape (3) Clamp (2) End Cap

IMPORTANT RISK INFORMATION

Indications For Use: The Medcomp[®] Titan HD is indicated for use in attaining Long-Term vascular access for Hemodialysis and Apheresis. It may be inserted percutaneously and is primarily placed in the internal jugular vein of an adult patient. Alternate insertion sites includes the subclavian vein as required.

Contraindications: This catheter is intended for Long-Term vascular access only and should not be used for any purpose other than indicated in these instructions.

Refer to Instructions for Use provided with the product for complete instructions, warnings, precautions, and contraindications. Observe all instructions for use prior to using products. Failure to do so may result in patient complications.